Smartline

Assistant 6000

Multi-Function Module for Smartline Systems

The **Assistant 6000** is a combination module, optionally fitted with up to three single modules (Pump 100 with or without pressure transducer, Smartline UV Detector 200 and/or valve drives with 2, 6, 12, or 16 switching positions).

The single instruments as well as the power supply are placed in one housing which can easily be integrated into a Smartline tower or a Prepline system.

Two of the instruments are situated at the front side whereas the third one, always an electrical valve drive, faces to the lefthand side.

The instrument operation and control is not changed by integrating into the Assistant 6000.

Counting only this instruments and not considering additional variations like the type of detector flow cell, the type of pump head or the valve specification (material or connectors) a total number of more than 300 variations is possible. This number itself already underlines the great variability of the Assistant 6000.

Thus, special modifications of the Assistant 6000 cover different segments of liquid chromatographic work, like

- ▶ Preparative Chromatography
- ▶ Solid phase extraction (SPE)
- ▶ Two dimensional Chromatography
- ▶ Column switching

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Modification Coding and Equipment Examples

Some of our system solutions are based on the Assistant 6000 for instance the standard preparative and SPE modules. – This list does not claim completeness. Within the above mentioned limitations we equip the Assistant 6000 totally according your selections.

Modification Coding of the Assistant 6000

The order No. of the Assistant 6000 is A5003 with the appendix Vxxx. The 1st digit describes the rear position, the 2nd and 3rd digits the front positions left and right. The single digits have the following meaning:

- 0 → Position not used
- 1 → Valve drive S6 with 2 switching positions
- 2 → Valve drive S6 with 6 switching positions
- 3 → Valve drive S12
- 4 → Valve drive S16
- 5 → Pump 100
- 6 → UV Detector 200
- 7 → Pump 100 with pressure transducer

Equipment examples of the Assistant 6000

In any modification of the Assistant 6000 the single instruments are included as basic modules. The desired pump head, flow cell, and valves are any way to order additionally to the Pump 100, UV Detector 200 and valve drives S6, S12, or S16 respectively.

Assistant 6000 V315 (Basic for Preparative Module 6100)

- 3 → Valve drive S12 with 12 switching positions (rear)
- 1 → Valve drive S6 with injection valve (left)
- 5 → Pump 100 without pressure transducer (right)

Assistant 6000 V415 (Basic for Semi-Preparative Module)

- 4 → Valve drive S16 with 16 switching positions (rear)
- 1 → Valve drive S6 with injection valve (left)
- 5 → Pump 100 without pressure transducer (right)

Assistant 6000 V022 (Selection Switching Module)

- 0 → The rear position remains unused
- 2 → Valve drive S6 with 6 switching positions
- 2 → Valve drive S6 with 6 switching positions

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Ordering Information

Order No.

A5003 V__*) **Smartline Assistant 6000, variant __*)**
Flow cells, pump heads, valves, and other accessories must be ordered separately.
*) three digit code of selected version

	Complete Systems
A5004 V03	Preparative Module 6100
A5004 V15	CHANCE HPLC System
A5004 V20	GPC Cleanup Unit 6500
A5004 V41	Sample Preparation Unit 6300
	Other complete systems on request

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Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.