
SmartMix

... more than just a mixer

SmartMix – more than just a mixer

The KNAUER SmartMix represents a new generation of static mixing devices, providing **outstanding mixing performance by applying a completely new type of mixing technology***. The mixing efficiency of SmartMix static mixers outperforms other commercially available static mixers and is even comparable to that of high quality dynamic mixers, all while featuring a **low mixing volume** design for improved gradient performance and short delay times. SmartMix static mixer combines several functions into one unit:

- ▶ a low pressure gradient mixer
- ▶ a high pressure gradient mixer
- ▶ optional flow splitting
- ▶ inline solvent filter

SmartMix static mixers is applicable for a **wide range of eluents** such as water/acetonitrile, water/methanol, and buffer gradients. With cartridges for both analytical (0.5–2.0 ml/min) and micro LC flow rates (0.1–0.5 ml/min or 1–100 µl/min with flow splitter respectively), the SmartMix mixer always provides the lowest dead volume. Just exchange the internal mixing cartridge to convert from analytical (350 µl) to the micro LC (<20 µl) version.

Since SmartMix technology is independent of the liquid flow direction, the SmartMix static mixer can be used to improve the performance of both low pressure gradient (LPG) or high pressure gradient (HPG) systems by simply inverting the flow direction. In addition, the SmartMix can also be used to establish **micro LC** flow rates by connecting a restrictor for flow splitting. The SmartMix mixer also includes an inline **filter** unit for added functionality.

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Key SmartMix benefits


- ▶ highly efficient mixing
- ▶ improved gradient performance
- ▶ significantly reduced baseline noise
- ▶ higher sensitivity
- ▶ low mixing volume
- ▶ short delay times
- ▶ versatility

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Low pressure gradient mixing mode

The portioning principle applied in low pressure gradient systems (namely, eluent “plugs” being formed in series) poses a **challenge** to the subsequent mixing device. SmartMix technology was designed to handle such eluent plugs efficiently to provide optimal mixing performance.

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
Highly efficient LPG mixing

The low pressure step gradients depicted here demonstrate the efficiency of different mixers for a demanding water/acetonitrile (ACN) eluent mixture.

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
High pressure gradient operation

By reversing the flow direction through the SmartMix static mixer, the mixing efficiency of high pressure gradient systems can also be easily **enhanced**.

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
Split mode operation for capillary LC

Use the micro LC cartridge of the SmartMix in conjunction with a restrictor (capillary) to achieve flow rates from 100 $\mu\text{l}/\text{min}$ **down to 1 $\mu\text{l}/\text{min}$** .

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

Less noise – higher sensitivity


The comparison below shows the **gain in sensitivity** by using the SmartMix static mixer for a reversed-phase (RP) separation carried out with a $\text{H}_2\text{O}/\text{ACN}$ gradient.

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linear gradient: $\text{H}_2\text{O}/\text{ACN}$, flow rate: 1 ml/min


High quality buffer/salt mixing

The buffer mixing example shown below demonstrates the **gain in baseline quality** achieved with use of the KNAUER SmartMix mixer in hydrophobic interaction (HIC) and other biochromatography applications.

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The performance of HPLC gradient systems can be significantly improved with the new SmartMix multifunctional static mixer. Thanks to its compact size and versatile design, the multifunctional SmartMix unit can also be easily integrated into almost any existing HPLC system.

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Technical Data

Flow range analytical: 0.5–2.0 ml/min;
 micro LC: 0.1–0.5 ml/min, 1–100 µl/min with flow splitting
Mixing volume analytical: 350 µl; micro: <20 µl
Pressure rating tested up to 400 bar (5800 psi)
Dimensions 29x45x29 mm (WxHxD)

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Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.

KNAUER · ASI · Advanced Scientific Instruments
 Wissenschaftliche Gerätebau
 Dr. Ing. Herbert Knauer GmbH
 Hegauer Weg 38 · D-14163 Berlin, Germany
 Telephone: +49 (0)30 / 80 97 27-0
 Telefax: +49 (0)30 / 8 01 50 10
 E-Mail: info@knauer.net
 Internet: www.knauer.net

Your local distributor: