**WellChrom K-2301/K-2401
Refractive Index Detectors
Refraktionsindexdetektoren
Manual / Handbuch**

V7149, 06/2007

CONTENT

Using this Manual	4
Conventions used in this manual	4
SOP's in this manual	4
Features of the WellChrom RI Detectors	5
Measurement Principle and Light Path	5
Calculation of the Signal Value	6
Installation of the RI Detector	6
Unpacking	6
Packing List.....	7
Protecting Yourself and Your Instrument.....	7
Front and Rear View of the RI Detectors.....	7
Function of Foil Keys	8
Grounding the Signal Cable.....	8
Connecting the KNAUER Interface Box	9
Connecting a Personal Computer and Digital Data Input.....	9
Use of the RI Detector Terminal Connectors.....	9
Assembling plug strips.....	10
Functions of the Control Connectors	10
Capillary Connection with a HPLC System	10
Operating the RI Detector.....	12
Power supply, ON/OFF, Self-Test	12
Internal Software Structure	12
The SIGNAL-Menu	13
The UNITS-Menu.....	14
The INT-Menu.....	14
The SIG Menu.....	14
The Menu $I_1 + I_2$	15
AUTOZERO and Zeroplate adjustment.....	15
The Menu CALIBRT.....	15
The FLUSH time Menu	16
The ANALOG OUT Menu	16
The GLP Menu.....	16
Error messages and their reasons	16
Cleaning and Maintenance.....	17
Terminate Operation and Shut Down	17
Removing Measurement Cell.....	17
Installing the Measurement Cell.....	19
Checking the Calibration.....	21
Spare Parts and Accessories	21
Technical Data.....	21
Warranty statement	41
Declaration of conformity.....	42

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Using this Manual

This manual refers to the WellChrom RI Detectors K-2301 and K-2401.

Conventions used in this manual

For easier reading, the term RI-detector is used for both instruments: the analytical detector K-2301 and the preparative detector K-2401. Only for specific information is the model number is mentioned.

Important Hints are marked by the marginal hand symbol.

Special Warnings are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates helpful advice.

SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks in the operation of your RI detectors. They include step-by-step instructions leading the user through all routine tasks. They can be used for documentation purposes and be copied, applied, signed, and filed in order to document the performance of the instrument.

Please operate the instrument and all accessories according to instructions and SOP's in this manual. This ensures accurate results and longevity of your equipment.

SOP 1	WAGO plug strip assembling	10
SOP 2	Emptying the measurement cell	17
SOP 3	Removing Measurement Cell of K-2301 RI Detector	17
SOP 4	Removing Measurement Cell of K-2301 RI Detector	18
SOP 5	Installing the Measurement Cell of K-2301 RI Detector	19
SOP 6	Installing the Measurement Cell of K-2401 RI Detector	20

Features of the WellChrom RI Detectors

The WellChrom K-2301 RI Detector is a high performance refractive index detector for use in analytical HPLC; the WellChrom K-2401 RI Detector is the corresponding instrument for use in preparative HPLC.

The instruments are controlled either by the installed firmware or externally by an interface box, connected to the control inputs. The device is suitable for laboratories functioning according to the FDA GLP regulations.

Measurement Principle and Light Path

A light beam emitted from a LED crosses the sample and the reference cell of the RI Detector twice, see Fig. 1. When both cells contain pure solvent, the system is calibrated to zero by means of a parallel zeroplate that positions the beam on the two detector diodes 1 and 2 in such a manner that the measured light intensities (I_1 and I_2) of the two diodes are virtually identical. When the sample cell contains a solution with a different refractive index, the light beam is geometrically proportionally deflected depending on the relative change of the refractive index (according to Snellius' refraction law).

This results in a change of the light intensities I_1 and I_2 , proportional to concentration and refractive index of the sample solution. From these intensity changes the signal value is calculated (as shown in the formulas in section „Calculation of the Signal Value“ on page 6) and indicated in the display.


Fig. 1 Light Path in the RI Detectors

The built in measurement cell of the K-2301 RI Detector operates with a measuring angle of 45° , the cell of the K-2401 RI Detector with an angle of 15° .

The 45° and 15° measurement cells are interchangeable with the effect that the sensitivity increases approximately by a factor of 3 using the 45° cell. The maximum flow rate will not be influenced by the measuring cell. It depends only on the inner diameter of the capillaries in the device (1.0 mm at K-2401, 0.3 mm or 0.7 mm at the K-2301).

The sensitivity and noise of the RI detectors are listed in the following table. For more technical details see page 21.

Table 1 Sensitivity of the RI Detectors

	K-2301	K-2401
Sensitivity	8×10^{-8} RIU	8×10^{-7} RIU
Noise	$\leq \pm 4 \times 10^{-8}$ RIU	$\leq \pm 4 \times 10^{-7}$ RIU

The detectors have a very low base line drift in the usable measurement range: $\pm 1 \times 10^{-3}$ RIU (K-2301) and $\pm 2 \times 10^{-3}$ RIU (K-2401).

Measurements are performed at a wavelength of $\lambda = 950 \pm 30$ nm. The signal detection and signal processing method allows the refractive index to be given out in the „online“ mode without absorption components. The display of the main menu also allows a permanent control of the detected light quantity. You will find additional information in section „Calculation of the Signal Value“. The autozero-range includes the complete measurement range.

By means of the RS-232 interface of the device, all measured data of the refractometer can be read out digitally. The KNAUER interface box enables use of analog control signals as well as analog data transfer. By using the KNAUER software packages ChromGate or EuroChrom, the control of the RI Detectors can be performed totally automatically.

Calculation of the Signal Value

The light beam reaches the two detector diodes 1 and 2 (see Fig. 1 on page 5) which deliver the intensity values during measurement, depending on the light beam's deflection. The results of the difference $I_1 - I_2$ and the sum $I_1 + I_2$ are calculated continuously and indicated in the menu SIGNAL.

The value of SIGNAL is calculated with the following formula:

$$\text{SIGNAL} = \left\{ \frac{(I_1 - I_2)}{(I_1 + I_2)} + a \right\} \cdot c \quad (1)$$

With $I_1 - I_2$ difference of intensity values
 $I_1 + I_2$ sum of intensity values
 a constant factor, defined by the function AUTOZERO
 c constant factor, defined by the function CALIBRT

The signal result, indicated in the menu SIGNAL, is given out to the device outputs.

The refractive index is highly temperature dependent. The change per 1°K for pure water is $\sim 1 \times 10^{-4}$ RIU and for typical organic solvents is $\sim 5 \times 10^{-4}$ RIU.

Installation of the RI Detector

Unpacking

After unpacking, please check the device and accessories thoroughly for any damage that may have occurred during transport. If necessary, put forward any claim for damages to the carrier.

Use the standard delivery list and check that the RI Detector delivery is complete. Please contact our service department if you are missing something or if you need support.

Packing List

The complete delivery of the RI Detector includes the following components:

- 1 WellChrom RI Detector K-2301 or K-2401
- 1 Euro power supply cable
- 1 Capillary connector with syringe
- 1 Operation manual
- 1 Plug strips for remote terminal control
- 1 RS-232 cable
- 1 coaxial connecting cable

All article numbers are listed in section Spare Parts and Accessories on page 21.

Protecting Yourself and Your Instrument

The use of the KNAUER K-2301 RI Detector as determined by the manufacturer, includes its usage in analytical HPLC systems. The use of the KNAUER K-2401 RI Detector includes its usage in preparative HPLC systems. The RI Detectors must be positioned at a vibration-free location. It should be placed away from draughts and intense heat, including direct sunlight.

For your own safety, you should observe the relevant regulations when working in a laboratory. Make sure that the main power supply is properly grounded and a proper three-pronged power cable is used. Connecting the RI Detector to faulty power sockets can cause injuries.

Front and Rear View of the RI Detectors

The front panel of the RI Detector is shown in Fig. 2. The rear panel in Fig. 3 shows all operating controls and connections.


Fig. 2 Front view of the RI Detector


The standard factory settings enable you to start the measurements immediately without the necessity of additional adjustments.

The RI Detector is completely controlled by the internal firmware. It can be operated by means of the keys on the front side.

Function of Foil Keys

The foil key area in Fig. 2 „Front panel of the RI Detector“ consists of four arrow keys and an AUTOZERO key.

Arrow Keys

Use the yellow arrow keys \triangleright „right“ or \triangleleft „left“ to move and to position the cursor on the display as well as to confirm entered or selected values.

Using the yellow arrow keys \blacktriangle „up“ or \blacktriangledown „down“ you can change the selected parameter or the available options of it.

AUTOZERO

The key AUTOZERO has four different functions:

- Pressing briefly in the menu SIGNAL will activate the „AUTOZERO function“, which sets the measurement signal and the integrator output to zero. This is performed by an appropriate selection of a new constant factor „a“ according to formula (1) on page 6.
- Extended pressing (> 3 s) in menu SIGNAL will activate a repositioning of the zeroplate (Fig. 1 on page 5). During this operation the message zeroplate moving is indicated.
- If the difference $I_1 - I_2$ exceeds a pre-set limit after power-on of the instrument, an automatic adjustment of the zeroplate is performed.
- Pressing this key at any site of the software causes an immediate return to the menu SIGNAL .


Fig. 3 Rear view and connectors of the RI Detector

Grounding the Signal Cable

The signal output of the RI Detector can either be grounded via the two outputs „GROUND“ or the shielding of the analog output. The device's ground is not connected to the signal ground. Signal grounding should be performed preferably via the analog output to the grounded integrator.

A multiple connection to ground (simultaneous grounding via „GROUND“ and analog output) has to be avoided. It could induce noise loops, that may lead to disturbed measurement results.

Connecting the KNAUER Interface Box

The KNAUER interface box is connected to a PC using its RS-232 interfaces. Data transfer from the detector to the interface box is enabled by connecting the analog integrator output of the RI Detector with the corresponding input on the interface box. For analog remote control of the RI Detector the box is connected according to the scheme in Fig. 4 and the information given in Table 2.


Fig. 4 Connections between the RI Detector and the KNAUER interface box

Table 2 Connections between the RI Detector and the KNAUER interface box

Socket at RI Detector	Socket at Interface- Box
INTEGRATOR 1V	INTEGRATOR INPUTS, 1
AUTOZERO	D.OUT3
FLUSH	D.OUT1
ERROR OUT	D.IN1
GROUND	D.GND

Connecting a Personal Computer and Digital Data Input

The RI Detector can be connected to the serial RS-232 interface of a PC using the RS-232 data cable. By means of an appropriate software, the data of the RI Detector can be digitally read out by the PC.

Use of the RI Detector Terminal Connectors

The terminal strip allows communication with further devices via electrical signals, also without using the interface box. Protect these connectors from electrostatic discharges to avoid damage of the electronics. Use suitable grounding equipment, or discharge yourself carefully prior to operating the connectors.

For wiring the Spectro Photometer K-2501 with other instruments you need cables with WAGO plug strips. Plug strips with 2, 3 or 4 positions are included in the standard delivery.

Assembling plug strips

They are mounted as follows:.


Fig. 5 Assembling plug strips

SOP 1 WAGO plug strip assembling

1. Insert the rounded end of the lever latch into the square opening of the selected connector of the plug strip.
2. Press the catch down as indicated by arrow.
3. Insert the uninsulated end of the cable into the opening under the catch.
4. Release the catch and remove the lever latch from the plug.

The cable is now firmly anchored in the plug strip.

Functions of the Control Connectors

The following functions are started while a signal >250 ms is applied to the corresponding input:

- AUTOZERO** A short circuit to GROUND starts the function Autozero. Measurements can be restarted only after completion of this process, indicated by the reappearance of the cursor on the display. Further information regarding AUTOZERO is given in section „Autozero and Zeroplate adjustment“ on page 15.
- FLUSH** A short circuit to GROUND starts flushing the reference cell with the pre-set duration of time. A permanent short circuit, duration >2 s, activates this function as long as the signal will be present.

The following signals are sent out by the RI Detector:

- INTEGRATOR 1 V** Integrator output, max. 1V bipolar with a resolution of 20bit
- ERR OUT** Error signal (TTL, open collector).

Capillary Connection with a HPLC System

The HPLC capillaries are connected to the fittings on the front side of the RI Detector, see Fig. 2 on page 7. The upper one is the input, the lower one the output.

Make sure that the capillary input and output connectors will not be confused; otherwise the measurement cell can be damaged.

The flow schemes of both RI detectors are shown in the following figures.

Connect the output of the HPLC column to the input of the RI Detector. Use a capillary with a low inner diameter (e.g. 0.3 mm) to keep the dead volume small.

Connect the capillary output of the RI Detector by means of a stainless steel or PTFE-capillary with a minimum diameter of 1.0 mm to the waste container.


Fig. 6 Flow scheme of RI Detector K-2301


Fig. 7 Flow scheme of RI Detector K-2401

Make sure that no back pressure is built up on the outlet side. This also may lead to damage of the measurement cell.

The maximum admissible flow rate through the K-2301 RI Detector is 5 ml/min

The maximum admissible flow rate through the K-2401 RI Detector is 100 ml/min

For more viscous solvents (e.g. 2-propanol) it is not recommended to work at the maximum flow rates.

Before running the solvent through the device make sure that the current solvent is completely miscible to the last one used. If not, the cell must be rinsed in an intermediate step with a solvent miscible with both ones.

Operating the RI Detector

Power supply, ON/OFF, Self-Test

The RI Detector operates with 115 or 230 V AC at 47-63 Hz. The setting is made by the manufacturer on customer's request. The standard setting is 230 V. If for any reasons it needs to be changed, please contact our service department. The actual setting is indicated on the rear panel of the instrument (Fig. 8).

Fig. 8 Indication of the input voltage

Make certain that the correct voltage has been set on the rear panel of the instrument, the power supply is grounded and a corresponding 3 prong power cable is used.

Connect the RI Detector to the power supply and switch it on with the power switch. Now the software release number appears in the display of the device:

Then a self test is carried out with checks of the following functions:

- Check of the integrator output for correct function and performance of an automatic offset-correction.
- Automatic adjustment of the zero plate.
- Automatic zero correction of the output signal (AUTOZERO).

Pay attention to the notes given on the display during self-tests. In case of error messages you can use the section Error messages and their reasons on page 16.

Prior to measurements, rinse measurement cell and reference cell with solvent and let the HPLC system get balanced to equilibrium, waiting a time period of at least 15 minutes. Then you can start measurements with your RI Detector.

Internal Software Structure

The software is divided in various menus, each of which allows particular settings and operational modes. You reach the single menus by positioning the cursor on the rhombus field \blacklozenge followed by pressing the \blacktriangle or \blacktriangledown keys. As indicated in Fig. 9, the menus will be called up in an endless loop.



Fig. 9 Menu sequence of the RI Detector

Now the single menus will be described in detail. Inside of any menu the cursor can be moved to the next or previous field using the \triangleright or \triangleleft keys respectively. There it is possible to increase or decrease the corresponding parameter settings by help of the Δ or ∇ keys. In some cases you can scroll with the Δ or ∇ keys through the available options. Moving to an other entry field using the \triangleright or \triangleleft keys respectively the entered changes will be confirmed.

The SIGNAL-Menu

The signal menu represents the main menu of the RI Detector and displays the following information:

- The measured signal intensity (signal)
- The measured light intensity is displayed in % if it falls below the pre-set value for **limit** in the $I_1 + I_2$ menu, otherwise this field remains free
- The set time constant in seconds

Additionally you can control in this menu the following functions:

Selecting the time constant

Using the time constant **t**: you can achieve a signal smoothing. It's value can be set to **0.1, 0.2, 0.5, 1, 2, 5** or **10** seconds. The larger this value is set the more the signal will be smoothed. A discrete change of input data will be recognized as full-scale deflection in the measuring signal after a duration of three times the time constant **t**. A time constant of 1 s fits most analytical purposes.

Flushing the Reference Cell

You can rinse the reference cell by placing the cursor to on "Flush" using the keys \triangleright or \triangleleft , and activate this function with the keys Δ or ∇ . The duration for flushing can be set in menu FLUSH, selecting item **time**. The standard setting is 30sec. Flushing is finished by pressing any key. An activated flushing process is indicated by the disappearance of the cursor in the display.

During flushing, measurement and reference cell are switched in series such that both cells are rinsed with solvent.

The flush flow rate should not exceed 20 ml/min. The flush solenoid valve and/or the reference cell may be destroyed when flushed with higher flow rates.

Signal Inversion

Effective for all outputs the sign of the signal value can be inverted when the cursor is placed on the menu item SIGNAL. Move the cursor with the keys \triangleright or \triangleleft and change the signal sign with the keys Δ or ∇ .

For the K-2301 RI Detector **SIGNAL** indicates a positive value and **signal** indicates a negative signal value. For the K-2401 RI Detector the situation is opposite: **signal** indicates a positive signal value while **SIGNAL** indicates a negative signal value.

The UNITS-Menu

UNITS	t:	s
◆	signal:	mRIU

This page gives information on the used units, Refraction Index Units (more precisely δ RIU), and units of time.

The INT-Menu

INT	I_1 :	0.4943
◆ 1/16	I_2 :	0.4881

In this menu page the measured intensity values I_1 and I_2 are indicated separately. Please note that the sum cannot be calculated directly from the indicated data, since these single intensities still contain an offset.

A continuously low intensity value can indicate an air bubble or a contamination in the cell. The cells can be rinsed using the function „Flush“.

The SIG Menu

SIG	$I_1 + I_2$:	0.9021
◆ ZP	$I_1 - I_2$:	0.0767

In this menu the values $I_1 + I_2$ and $I_1 - I_2$, calculated from I_1 and I_2 , are displayed. Best measurement results are achieved if both I_1 and I_2 are close to 0.5, which will result in $I_1 + I_2 \sim 1.0$ and $I_1 - I_2 \sim 0$.

Furthermore, a zeroplate adjustment can be activated by selecting the function **ZP**.

The Menu $I_1 + I_2$

$I_1 + I_2$: on line ◆ limit: 80 %
--

In menu $I_1 + I_2$ the following items can be selected:

- online:** $I_1 + I_2$ is continuously calculated for every current data point. Thus it is ensured that fluctuating absorption components have no influence on the measurement value. The noise component in the measurement signal increases slightly.
- AUTOZERO:** $I_1 + I_2$ will be kept constant during the entire measurement. The current intensity of the sum is recalculated after activation of AUTOZERO and saved.
- fixed ->1:** For the entire duration of measurement $I_1 + I_2$ is defined as 1. This mode corresponds to the detection method of the former KNAUER RI 9800 types.
- limit:** By selecting the item limit, the use of the keys Δ or ∇ allows a percentage to be set. If the detected light quantity falls below this limit, it will be displayed in the menu Signal. Otherwise no light quantity data are displayed.

AUTOZERO and Zeroplate adjustment

The best measurement results are obtained, if the sum of $I_1 + I_2$ is close to 1 (corresponding to a detected light intensity of 100%) and if the difference $I_1 - I_2$ is close to zero. The minimum detected light intensity ($I_1 + I_2$) should amount to more than 80%. If it reaches less than 80%, the baseline will get more noisy. The intensity values I_1 and I_2 are displayed in the menu **INT**.

The electronic function Autozero serves for correction of a base line offset due to drift, that may be caused by thermal or solvent fluctuations.

Large deviations from the optimal values can be corrected by a zeroplate adjustment. First it should be checked if deviations in I_1 and I_2 are being caused by an air bubble in the cell. This can be recognized by a very small light intensity $I_1 + I_2$ (Display in menu Signal) or by a big difference $I_2 < I_1$. An air bubble in the cell can be removed by rinsing with solvent and by using the Flush function to flush the reference cell.

During an error-free operation the zeroplate adjustment can be activated by selecting the page **ZP** in the menu **SIG** or by pressing the key Autozero for more than 3 s.

The Menu CALIBRT

In this menu the calibration constant c according to formula (1) on page 6, to be used for signal value calculation can be selected. Available settings are in the range $0.125 \dots 8.000 \times 10^{-3}$ for instruments with software releases up to V1.4, and the range $0.125 \dots 128.0 \times 10^{-3}$ for instruments with software release V1.5 and higher. The current software release can be read on the display after switching on the instrument.

For changing the setting place the cursor on the position **coarse** or **fine** using the keys \triangleright or \triangleleft . Now you can select the desired value by means of the Δ or ∇ keys. With **coarse** the values are doubled or divided by 2, with **fine** you can change the digits in steps of 4 at the third decimal position.

All KNAUER detectors are adjusted to a standardized sensitivity in-factory by comparing measurements with calibrating solutions at $T=20^\circ \text{C}$. The selection of another value for calibration applies to all measurement outputs, and erases/overwrites the previous entry.

Changing the calibration constant c does not influence the signal-to noise ratio! By selecting this factor you can achieve any calibration for the RI Detector. Thus, an absolute calibration by means of a solution with known concentration at a certain temperature can also be performed.

The FLUSH time Menu

FLUSH time: ◆ 030 s

You can define in this menu the time duration for flushing the reference cell. The available range is 1...900 s. Place the cursor with the keys \triangleright or \triangleleft to the time field, and select the desired duration with the keys Δ or ∇ . While pressing the ∇ key you can reach the function **no valve** beyond the value 0, which is used to deactivate the valve circuit.

Flushing periods of 10...20 s are recommended for flow rates of approx. 1ml/min in the RI Detector. With lower flow rates the use of longer flushing periods is advisable.

The flush flow rate should not exceed 20 ml/min. The flush solenoid valve and/or the reference cell may be destroyed when flushed with higher flow rates.

The ANALOG OUT Menu

In the menu ANALOG OUT the scaling for the analog output signal can be chosen. The signal voltage of 1 V can be adjusted in 16 steps from 10^{-5} to 10^0 RIU.

The GLP Menu

In the menu GLP you get the serial number (S/N) of the instrument and in the second line the overall operating time is indicated.

Error messages and their reasons

The indication of one of the following error messages refers to a malfunction of the device:

Error Moving Zero Plate This message can appear if the zero plate adjustment upon switching on or after manual activation has failed. The reason can be an air bubble in the measurement cell or some other malfunction. Make sure that the measurement cell is flushed with solvent. Check the intensity values I_1 and I_2 (indicates presence of air bubbles) and regard the notes in section „Autozero and Zero plate adjustment“ on page 15.

Should this error message persist even after eliminating simple error sources or after multiple switching off and on again, please contact our service department.

Signal Overflow Overflow of one of the two AD (Analog-digital)-converters. Flush measurement- and reference cells again and press the key AUTOZERO.

Cleaning and Maintenance

If the light intensity falls below a preset value this will be indicated in the main menu. The actual intensity will then be shown in %. If you are working in the **on line** mode (menu $l_1 + l_2$) the light intensity will not influence the measuring value. However, the noise will increase with decreasing light intensity.

A decrease of light intensity may have different causes:

Problem	Solution
IR-absorption of eluent	use other eluents
Air bubbles in the measurement cell	flush the cell
measurement cell dirty	clean / exchange the cell
connection leaking, cell outside wet	remove the cell, clean it outside, reinstall it properly with new sealings

Terminate Operation and Shut Down

To terminate the operation of the RI Detector for longer periods, the measurement and reference cells must be emptied according to the following SOP:

SOP 2 Emptying the measurement cell

1. Activate the Flush function in the menu Signal. Make sure that the selected duration for flushing is long enough to complete the steps 2 to 5. Activate the function flush repeatedly, if necessary.
2. Fill solvent into the measurement cell using a 10 ml syringe
3. Carefully press air into the measurement cell using a 10 ml syringe
4. Fill acetone into the measurement cell using a 10 ml syringe
5. Carefully press air into the measurement cell using a 10 ml syringe
6. Close the connections at the front panel of the instrument by means of suitable fittings.

Buffer solutions can damage the detector cells after only a few hours. Therefore make sure that no buffer solution will remain inside the cells after shutting down the RI Detector, even overnight.

Removing Measurement Cell

Removing the measurement cell is required for cleaning or for exchange. Removing and reinstalling the cell should be done by authorized service personnel only, although it can be performed by trained users.

Remove the power plug at the rear side of the detector before opening the housing! Inside the instrument electrical connections are located carrying high voltage which may damage your health.

SOP 3 Removing Measurement Cell of K-2301 RI Detector

1. Empty the measurement cell according to SOP 2.
2. Loosen the four housing screws on the side panels and the two screws on top of the housing using a cross tip screw driver.
3. Take off the housing pulling upwards.

4. Locate the measurement cell inside the housing. The measurement cell is installed on the optical bench from above and is fastened as shown in Fig. 10 and Fig. 12 by means of a cap nut.


Fig. 10 Site of measurement cell in the K-2301 RI Detector

5. Loosen the cap nut.
6. Pull out the cap nut together with lock washers and cell cover.
7. Pull the measurement unit out of the housing.

You can clean the measurement cell in a ultrasonic bath filled with a suitable cleaning solution, or install a new one.

SOP 4 Removing Measurement Cell of K-2401 RI Detector

1. Empty the measurement cell according to SOP 2.
2. Loosen the four housing screws on the side panels and the two screws on top of the housing using a cross tip screw driver.
3. Take off the housing towards the top.
4. Locate the measurement cell inside the housing. The measurement cell is installed on the optical bench from above and is fastened as shown in Fig. 11 and by means of four spring-loaded screws.


Fig. 11 Site of measurement cell in the K-2401 RI Detector

5. Loosen the four recessed head screws with springs with a cross tip screw-driver.
6. Pull out the four screws together with springs and cell cover.
7. Pull the measurement unit out of the housing.

You can clean the measurement cell in a ultrasonic bath filled with a suitable cleaning solution, or install a new one.

Installing the Measurement Cell

Make sure that all parts are properly functioning before installing.

It is recommended to use a new set of seals whenever installing a cell.

SOP 5 Installing the Measurement Cell of K-2301 RI Detector

1. Put the lower seals over the capillaries' ends in the cell ground, considering the bevels.
2. Insert the measurement cell, considering it's position (see Fig. 12) into the housing. The writing KNAUER has to be situated on the upper side.
3. Put the upper seals over the capillaries' ends in the cell cover, considering the bevels. Insert the cover, having the flat side directed toward the instrument's front. Put in the lock washers and screw the cap nut into the optical bench.
4. Screw tight the cap nut with the hand.

Tightening the cap nut too weakly can result in leakage of the measurement cell. Tightening the cap nut too tightly can lead to a direct contact of the capillaries to the cell, resulting in damage to the cell.


Fig. 12 Layout of 45° measurement cell (left); and orientation on the optical bench (right)

- Put the housing over the instrument, and tighten the four screws.

SOP 6 Installing the Measurement Cell of K-2401 RI Detector

- Put the lower seals over the capillaries' ends in the cell ground, considering the bevels.
- Insert the measurement cell, considering its position (see Fig. 13) into the housing. The writing KNAUER has to be situated on the upper side.
- Put the upper seals over the capillaries' ends in the cell cover, considering the bevels. Insert the cover, having the flat side directed toward the instrument's front.
- Put in the springs and the four screws, and screw them into the optical bench .


Fig. 13 Layout of 15° measurement cell (left); and orientation on the optical bench (right)

- Screw the four recessed head screws tightly.

Tightening the screws too weakly can result in leakage of the measurement cell. Tightening the screws too tightly can lead to a direct contact of the capillaries to the cell, resulting in damage to the cell.

6. Put the housing over the instrument, and tighten the four screws.

Checking the Calibration

The RI detectors are calibrated by the manufacturer using an aqueous solution of citric acid (10mg/ml) towards de-ionized water, resulting in a measured value of 1.147×10^{-3} RIU at a room temperature of about 23°C. An aqueous solution of glucose monohydrate (8.28mg/ml) can also be used for calibration. This results in a relative refractive index of 1×10^{-3} RIU at 20°C.

Spare Parts and Accessories

15°- Measurement Cell	A0287
45°- Measurement Cell	A0294
Set of seals	A0277
Power supply cable 230 V	M1642
Power supply cable 115 V	M1651
RS-232 data cable (female - female)	A0895
Connectors, plug strips, signal cable 12 lines	A1402, A1404
Cable for integrator (cinch-cinch)	A1589
Cable for integrator (cinch-banana-plugs)	G1023
Manual	V7149

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Warranty statement

The warranty period of the WellChrom RI-Detectors K-2301 und K-2401 is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Defective RI-Detectors should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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Declaration of conformity

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Manufacturer's name and address:
Herstellername und -adresse

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comply with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)
EN 6100-3-2 (10/2006)
EN 61326-1 (10/2006)

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Dr. Ing. Herbert KNAUER GmbH