User Manual / Handbuch

Fast scanning Spectro-Photometer K-2600

Schnell scannendes Spektralphotometer K-2600

V7016, 09/2002

CONTENTS

Using this Manual.....	5
Conventions in this manual	5
SOP's in this manual.....	5
General Description	6
Basic Concepts of the K-2600 UV detector.....	7
Optical Path of the K-2600 UV detector.....	7
Preparing the K-2600 UV detector for Operation	7
Unpacking.....	7
Power Supply, ON/OFF, Self Tests	7
Installation of the Flow Cell.....	8
Front Panel Controls	8
Front View of the K-2600 UV detector	8
Function of Foil Keys	9
Arrow Keys	9
Numeric Keys	9
AUTOZERO.....	9
SCAN	9
Rear Panel View of the K-2600 UV detector	9
Operating the Photometer	10
Display and Menu Structure.....	10
Main Menu	10
Appearance of Cursor and Data Input.....	10
The SETUP Menu	11
Das GLP Menu	14
Das LINK Menu	14
The VIEW Menu	15
Direct Control of the K-2600 UV detector	15
Selecting wavelength(s)	16
Programming the K-2600 UV detector.....	17
Creating Programs	17
Program Execution.....	19
Creating Program Links	20
Programming a link	20
Running program links	21
Deleting Programs and Program Links	21
Wake Up-Program.....	21
Using the Wake Up program	22
Scanning UV-Spectra.....	22
The Scan Menu	22
Scan Autozero	23
Performing wavelength scans.....	23
Output of Scan Data.....	23
Connecting other Instruments to the K-2600 UV detector.....	24
Using the Event and Remote Sockets	24
Assembling Plug Strips.....	24
Event Socket	24
Remote Socket.....	25

RS 232 Serial Interfaces	25
The KNAUER Net	26
List of display messages.....	26
Software Control of the K-2600 UV detector.....	26
Preparing for Software Control.....	27
Installing the KNAUER Net	27
Configuring the KNAUER Net.....	27
Simple Maintenance	28
Control of the lamp's functionality.....	28
Changing the lamp	28
Cleaning the flow cell.....	28
Analytical flow cells	29
Preparative flow cells	29
Adjusting the path length of the preparative flow cells	30
Flow cells with fiber optical connectors.....	30
Spare parts and accessories	31
Standard delivery.....	31
Spare Parts	31
Available Flow cells for your K-2600 UV detector	32
Analytical Cells	32
Preparative Cells	32
U-Z View™ Mikro-Durchflusszellen / CE Zellen	32
Spare parts for flow cells	32
Analytical Cells	32
Preparative Cells	32
Technical Data	33
Declaration of conformity	34
Guarantee statement	35

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Using this Manual

This manual refers to the WellChrom K-2600 UV-Detector Firmware Revision 1.06 or higher. It is valid for any combination with analytical flow cells order number A 4061, A 4042, A 4045, preparative flow cells order number A 4066, A 4067, A 4068, A 4069, A4095 and all U-Z View™ micro flow cells in standard and fiber optic version.

Conventions in this manual


Special Warnings are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates **helpful advice's**.

Important Hints are marked by the marginal hand symbol.

Arrows in an outlined form like this example used in block diagrams, indicate an automatic program run and change to the next line without the necessity of manual interventions.

Arrows like this, used in block diagrams, indicate that the user is asked to press the corresponding arrow keys

The triangles symbolize the use of corresponding arrow keys.

SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks in the operation of your WellChrom K-2600 UV detector. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and be copied, applied signed, and filed in order to document the performance of the instrument.

Please operate the instrument and all accessories according to instructions and SOP's in this manual. This ensures proper results and longevity of your equipment.

SOP's in diesem Handbuch

SOP 1	Installation of the Flow Cell.....	8
SOP 2	Time Programming of Absorption Wavelength.....	17
SOP 3	Programming events for control of a fraction collector.....	18
SOP 4	Assembling Plug Strips.....	24
SOP 5	Changing the lamp.....	28
SOP 6	Cleaning an analytical flow cell.....	29
SOP 7	Cleaning a preparative flow cell.....	29
SOP 8	Changing the path length.....	30

General Description

The WellChrom K-2600 UV detector is a new, programmable and fast scanning detector. It is capable of measuring up to four wavelengths simultaneously. In stand alone mode two wavelengths can be used simultaneously and four spectra can be stored in the detector. The scan function provides you with spectral information of your substances over the whole wavelength range without stopping the eluent flow of your HPLC system.

Up to 20 time programs with 99 program lines can be stored and are permanent available, up to nine links between programs can be defined. The instrument is ideally suited for use as a stand alone instrument as well as for use in a complete system under control of the KNAUER HPLC software EuroChrom for Windows.

Wavelengths between 190 and 740 nm are selected using a moveable grating monochromator with an accuracy of ± 1 nm in 1 nm steps. The bandwidth is 8 nm or better. Additionally, a cut-off filter can be activated. This filter suppresses secondary grating orders when measuring with a deuterium lamp at wavelengths > 380 nm and thus ensures the selectivity of the measurement. A halogen lamp for the visible wavelength range can easily be installed by the user.

Data acquisition and control of the K-2600 UV detector are performed digitally and such allow trouble free operation of the instrument. Two adjustable analog outputs and the option for external analog control are also available. Additional features offered by the K-2600 UV detector are the control of other instruments via event outputs, the control of a fraction collector by time window and/or level sensing and the possibility to give out ratio plots, e.g. for peak purity control, via the analog outputs.

An integrated self-test and self-calibration performed after switching on ensures reliable and reproducible operation. The highly sensitive detector features a extremely low noise level ($\leq 1 \times 10^{-5}$ AU) and baseline drift ($\leq 10 \times 10^{-5}$ AU/h). The autozero function is available across the entire measuring range. The K-2600 UV detector satisfies all GLP requirements of a regulated laboratory.

A rich palette of flow cells available for the KNAUER photometers, ranging from cells for Nano-HPLC with flow rates < 100 nl/min to preparative cells for up to 10 l/min, make the photometer K-2600 highly flexible in the full range of LC applications.

Moreover, a version of the K-2600 UV detector equipped with fiber optical connectors of the cell is available. This enables the spatial separation of detector and flow cell and thus the use of these instruments in hazardous locations.

Basic Concepts of the K-2600 UV detector

Optical Path of the K-2600 UV detector


Fig. 1 Optical path of the K-2600 UV detector

The light emitted from the lamp [1] is bundled and guided to a grating monochromator [5] for defining the desired wavelengths. After passing a filter [7] (optional) the beam is splitted by a semi transparent mirror [8]. One beam gives a reference signal [11], the other beam is guided through the cell [9] where optical absorption is measured [10]. The difference or the ratio of the two signals can be calculated.

Preparing the K-2600 UV detector for Operation

Unpacking

After unpacking, please check the device and accessories thoroughly for any damage that may have occurred during transport. If necessary, put forward any claim for damages to the carrier.

Use list "Standard delivery" on page 31 and check that the K-2600 UV detector delivery is complete. Please contact our service department if you miss something or if you need support. Please fill out the guarantee registration card and return it to us immediately.

Power Supply, ON/OFF, Self Tests

The K-2600 UV detector is equipped with an universal power supply which operates in the range from 90 to 260 Volts AC. A manual setting of the supply voltage is therefore not required.

Make sure that mains power supply is properly grounded and a correct three-pronged power cable is used.

Connect the power cable to the device and the socket, and then switch the device on with the ON/OFF key, see [5.2] and [5.3] in Fig. 3 on page 8.

After switching on, the display, see 1 in Fig. 2 on page 8, is indicating some information, such as **HPLC SPECTRAL PHOTOMETER**, and the **VERSION** of the internal Software, e.g. **1.07**.

The power-on selftest performed includes a check of the electronic function, wavelength calibration, see „AUTO CALIBRATION“ on page

11 and an autozero. Successful completion of selftest is indicated by **systemtest ok**.

Then the lamp of the device is heated up until the final, constant working temperature has been reached. During heating up the display shows **hea** in the lower left position. The starting routine is completed by an automatic calibrating procedure, and when indicating **on** the device is ready for operation.

Installation of the Flow Cell

The K–2600 UV detector comes factory configured without flow cell. Before operating the photometer, you need to install a appropriate flow cell from the palette of flow cells available from KNAUER.

SOP 1 Installation of the Flow Cell

This instruction refers to the K–2600 UV detector without fiber optical connectors. Refer to Fig. 2 „Front panel of the K–2600 UV detector“ on page 8.

1. Loosen the two „Knurled cell screws“ [3] with the fingers and remove them.
2. Pull out the „Flow cell housing“ [4].
3. Take the cell or dummy cell with two fingers and remove it upward.
4. Insert the new „Flow cell“ [7] and make sure that the engraved specifications point towards the user (...can be read), and the fixing hole on the back side of the cell meets the corresponding metal pin of the photometer’s housing.
5. Now push the complete system towards the housing, insert the two screws and thighten them manually.

Front Panel Controls

Front View of the K–2600 UV detector

1. Display information
2. Foil key area
3. Knurled cell screws
4. Flow cell housing
5. Flow cell cradle
6. Outlet
7. Flow cell
8. Path length of cell
9. Inlet


Fig. 2 Front panel view of the K–2600 UV detector

Function of Foil Keys

See „Foil key area“ [2] in Fig. 2 „Front panel of the K–2600 UV detector“.

Arrow Keys

The yellow, arrow-shaped keys can be used for moving and positioning the cursor in the display and for confirming an entered value.

Numeric Keys

All ten numeric keys 1 to 0 enable entering of numeric values at the position of the cursor.

In certain menus, the keys 0, 1 and 2 are used to RUN (1) or to HOLD (0) a defined program, or to insert event programmings, e.g. 2 to program the event pulse \square .

AUTOZERO

When pressed shortly, this performs an auto-zero over the whole wavelength range. Pressing longer than two seconds starts a scan autozero.

SCAN

The yellow colored SCAN key switches to the SCAN menu. After switching to this menu, pressing the SCAN button again triggers a scan over the chosen wavelength range.

Rear Panel View of the K–2600 UV detector


Fig. 3 Rear panel view of the K–2600 UV detector

5.1	Fuse compartment	5.6	Fan
5.2	Power connector	5.7	Serial number
5.3	ON/OFF switch	5.8	Terminal strip
			Remote connectors
5.4	RS-232 interfaces	5.9	Terminal strip
			Event connections
5.5	Socket for Ethernet-connection	5.10	Analog outputs
			(to recorder or integrator)

Operating the Photometer

Display and Menu Structure

After switching on the instrument, see chapter „Power Supply, ON/OFF, Self Tests“ on page 7, the device is ready for operation, showing on the display the main menu or operation screen. The center part of Fig. 4 shows an example of the main menu.

Using one of the arrow keys of the foil keypad, the blinking cursor can be moved to all positions of the menu. If the cursor has reached an exterior position in the menu, e.g. **ON** in the lower left corner of the main menu, one more pressing of \triangleleft for at least 2 seconds will indicate the VIEW MENU on the display.

Main Menu


Fig. 4 Overview menu structure

Some menus consist of several pages, indicated by a rhombus \blacklozenge in the lower left corner of the display. If you place the cursor on this rhombus, all menu pages can be accessed by pressing \triangle or ∇ .

Appearance of Cursor and Data Input


Fig. 5 Blinking cursor Cursor waiting for data input

The blinking cursor can be moved to any desired position, see Fig. 5 (left).

Once the position for data input has been reached, simply press the first numeric key to start data input. Now the cursor changes to underline character beside the first input, as you can see Fig. 5 (right). Enter next digit and proceed until numeric input is complete. Confirm and leave this position by pressing any arrow key.

Data input mode is disabled automatically, if a time limit of 60 seconds without pressing any key is reached.

The SETUP Menu

In the SETUP menu parameters for the operation of the K-2600 UV detector can be adjusted. These parameters comprise values for operation of the instrument.


Fig. 6 Block diagram overview SETUP menu

TIME CONSTANT

You can select smoothing the signal with a time filter. The value can be chosen between **0,1; 0,2; 0,5; 1; 2; 5** and **10** seconds. The bigger the value for the time constant, the more the signal is smoothed. For most analytical purposes, a time constant of **1s** is best suited.

AUTO CALIBRATION

(Second line in the window TIME CONSTANT):

Moving to this field and pressing or starts the procedure for the calibration of the wavelength scale. It searches automatically for reference points within the spectra of the lamp for reliability and reproducibility of the wavelength. The points used for calibration are the zero order reflection of the grating (000 nm) and the H α spectral line of Deuterium (656 nm). Auto calibration is performed automatically on power on.

LAMP

You can specify the lamp type (**Deuterium/Halogen**). Changing the lamp type resets the lamps working time after a control question, the lamp number is incremented. The number of the lamp cannot be edited. In the field with working time of the actual lamp you can reset this GLP value after a control question. The lamp number is then incremented. This function should be used only when replacing a deuterium lamp. The **ON/OFF** field enables you to switch off the lamp without switching off the instrument.


Fig. 7 Menu for lamp specifications

Signal

You can select the signal source used for all data outputs. This includes display on the operation screen as well as analog and digital data outputs. You can select the options **absorption**, **signal** and **reference** channel. The values **sig** and **ref** give the measured light intensity on the signal and the reference channel respectively. The sig and ref values range between 0 and 1, their absolute values are dependent on the wavelength, the solvent, the flow cell used and the age of the lamp.

SCAN OUT

You can select the output speed of scan data via integrator output. Available values are **1, 2, 5, 10, 20** and **50** steps per second.

ANALOG OUT FULL SCALE

You can scale the two integrator outputs according to your needs.


Fig. 8 Menu for scaling the analog outputs

In the first field you can choose the signal range covered by the output signal. Values that can be chosen by arrow keys from 10^{-4} to 10 AU in 16 steps. In the right hand field you can choose the maximum output voltage from **0.1, 1** and **10 V**.

ANALOG OUT OFFSETS

You can set an offset to the two integrator outputs. The output voltage of the outputs is set to zero during the self test of the instrument. In this window, you can set a numerical value for an offset in mV units that will be added to the corresponding signals on the analog outputs. This function can be useful when using a chart recorder to separate the charts for the two channels.

ANALOG IN

You can calibrate the analog input for external wavelength control.

Activating Set-Zero by pressing Δ or ∇ on this field defines the actual input voltage as zero, wavelength 000 nm. Then a signal is applied, e.g. 5V. The actual „λ“ field shows the interpreted wavelength, that should

be 500 nm with the preset conditions of 100 nm=1 V. This wavelength value can then be modified with the arrow keys Δ or ∇ , which creates a second calibration point.

CONTROL

The instrument can be set to different external control modes. You can choose between **KNAUER-NET**, where the baud rate is currently limited to 9600, **RS-232**, **RS-485** and **ANALOG**. When choosing the **ANALOG** option, only one wavelength can be used and controlled by the analog input voltage. The operation screen is modified correspondingly. The options RS-232 and RS-485 enable the use of a RS-232 or a RS-485 protocol for controlling the instrument with a terminal program.

The option **KNAUER-NET** is used when operating the K-2600 UV detector as a stand alone instrument or for combinations with other KNAUER instruments, e.g. for external control with the EuroChrom software (see chapter „Software Control of the K-2600 UV detector“ on page 26.

DATE/TIME

Allows setting system date and time.

FRACTION

This window provides you with the possibility to control a fraction collector by using the event outputs of the K-2600 UV detector. Additionally, signal levels to be exceeded for fraction collection and a delay time respecting the dead volume of the capillaries between detector and collector can be set.


Fig. 9 Menu for setting up the fraction collector control

Enter a delay time in seconds in the first line. Set the events for control of the in the second line, the event that controls the fraction collector can be chosen. Possible values are **E1** or **E2**.

E1 activates a relay output, E2 a TTL output. The corresponding connector on the rear side of the instrument can be connected to the fraction collector.

In the field **level** you can set a minimum signal level that must be exceeded before fraction collection starts. Both functions, **level** and **event**, are evaluated in a logical **and** functions. This means, that the E1 output is activated when the level is exceeded **during** activation in the time table of a method. In case of missing time programs the appropriate event has to be set to „1“ in the operation screen. If fraction collection is controlled by the signal level only, the event output becomes active every time the signal exceeds the preset signal level.

Das GLP Menu

The GLP report reports the operating statistics of the photometer.


Fig. 10 Block diagram GLP menu

Das LINK Menu

The link menu is used when running links, i.e. the pre-defined run of program combinations. In this menu the momentary status of the program links can be monitored.


Fig. 11 Block diagram LINK menu

In the first line of the menu the momentary status can be observed, the number of runs **R-xx** to be performed counts downwards. In the second line the next step of the link to be executed is shown. Arrows at the right side of the display indicate, that more program lines above and below, respectively, do exist.

For more information on links, see chapter “Creating Program Links” on page 20.

The VIEW Menu

The View menu gives an overview of existing programs and program links. See Fig. 12 for an example where programs 01, 02, 03, 04 and 11 and link 21 have been defined and do exist.


Fig. 12 Block diagram VIEW menu

Direct Control of the K-2600 UV detector

Direct Control of the K-2600 UV detector is performed via operation screen. Many functions like choice of wavelength(s) except scanning, can be controlled from the operation screen.

Programs can be used for determining and changing wavelengths or setting signal events during runs. The K-2600 UV detector features storage and selection between up to twenty different programs with a maximum number of 100 program lines for each program. You can store nine different links for calling programs within programs.

The most powerful feature of the K-2600 UV detector is its ability to scan UV-spectra with a speed of approx. 100 nm/s on flow. Scans are handled in the SCAN menu, which is accessed simply by pressing the SCAN button.

After power on, the instrument starts automatically with the last used program. In case no programs had been loaded, the standard program is defined to be No. 01, as shown in Fig. 13.


Fig. 13 Standard display with cut-off filter activated

Selecting wavelength(s)

Position the cursor on one of the two wavelength fields and select the wavelength by entering the desired values with the numeric keys. If you want to use one single wavelength, enter a Zero in the unused fields. Entering zero deactivates the respective wavelength. If you use both wavelengths, the polychromator moves between the positions corresponding to the selected wavelengths. This movement is performed with a frequency of 10 Hz that cannot be changed and. You can hear the motor of the polychromator as a high frequency tone.

A **cut-off filter** for suppressing secondary grating orders is activated automatically if one or both wavelength(s) exceed 381 nm and a deuterium lamp is installed. The use of the filter is displayed by a superscripted ⁺, as shown in Fig. 13.

Output options: The Output of absorption signals can be configured. Move the cursor to the signal fields. Choose between different options by pressing any numeric key. Options are:

- Signal inversion indicated by a little superimposed minus sign.

Two different ratio plots are available:

- Wavelength (1) divided through wavelength (2) and
- vice versa.

Ratio plots are indicated by symbol 1/2 and respectively 2/1 at the decimal point of the absorption value Fig. 14.


Fig. 14 Inversion of the signal and ratio plots: Set indicated output options

Ratio plots contain information on peak purity. It is calculated by dividing the intensity values at two different wavelengths. It is possible to measure the absorption at one wavelength and to monitor simultaneously the ratio of absorption between this wavelength and another one. Ratio plots can be accessed at either analog out connector or via digital data acquisition.

The scaling of the analog outputs can be adjusted to different values of **Absorption Units** or to different voltage values. For details see section „ANALOG OUT FULL SCALE:“ on page 12 in section „The SETUP Menu“.

EVENTS — The EVENT outputs on the rear side of the instrument provide time programmed signals to other instruments. All features, together with that of the remote control socket, is described in chapter „Connecting other Instruments to the K-2600 UV detector“ on page 24.

When working without a time program, using events enables fraction collection, as described in section "FRACTION" on page 13 in chapter „The SETUP Menu“.

Programming the K-2600 UV detector

Creating Programs

Time displays are given in minutes with decimal calculation of seconds, i. e. 0.3 min = 18 s. Time resolution during input is 1/10 min, during program runs 1/100 min. All programs are saved permanently.

Wavelengths and events can be changed as a function of time. Setup parameters such as scaling of the analog outputs cannot be time programmed. All options for analog outputs like ratio plots can be defined for entire programs and will be used when loading this program.

The K-2600 UV detector can store up to twenty programs. By moving the cursor to the VIEW menu, you can check the availability and allocation of program numbers.

Use the EDIT mode for modifying or creating programs. Move the cursor to the „**Edit Prog**” field and enter the desired program number with numeric keys. Any numeric values can be overwritten. Pressing any arrow key confirms the actual value.

The first line of a program is the start time at 0 min. This time cannot be changed. Message **Start time is fixed** appears on editing this line. The maximum value for time is 999 min.

To edit a new time line, move the cursor to the time field and press ▽. If this line is empty, in the time field will be displayed *****.***


Fig. 15 Empty line in program no. 3

You can delete a program line if you enter 0 in the field time and press any arrow key. The message **Delete this line? Confirm by cursor** will be displayed. You can confirm deletion if you press any arrow key. If you press a numeric key, or wait for 60 seconds, deletion will be aborted.

Values 0, 1 and 2 can be used for programming of events, see chapter „Connecting other Instruments to the K-2600 UV detector“ on page 24.

SOP 2 Time Programming of Absorption Wavelength

Use this instruction for programming the absorption wavelength as a function of time.

1. Activate the operation screen, see Fig. 4 on page 10.
2. Move cursor to the field Edit Prog, press any numeric key to generate a numeric label for the new program. Press any arrow key to confirm and enter edit mode for the new program.
3. The cursor is automatically moved to Time 000.0, which is the fixed start time.

4. Press \triangleright to go to field $\lambda 1$. Enter desired wavelength. Move cursor to field $\lambda 2$ and enter your second desired wavelength. The display should resemble section B in Fig. 16 on page 19.
5. Press ∇ to create a new program. The display should resemble section C in Fig. 16 with the blinking cursor on the first asterisk.
6. Press 1 and 0 and then \triangleright to define time step 10.0 min. Enter field $\lambda 1$ and put in 240 nm and complete entering by pressing ∇ . The value for $\lambda 2$ (254 nm) is kept constant as long as no changes are made.
7. Repeat the last two steps, define time step 20.0 min and enter 280 nm. The display should resemble section D in Fig. 16.

To run a program it must be loaded, see chapter „Program Execution“ on page 19.

SOP 3 Programming events for control of a fraction collector

Use this instruction to control a fraction collector via the K-2600 UV detectors EVENT socket using signal level **and** a defined time window.

1. Enter the SETUP menu and proceed to the window FRACTION, see Fig. 9 on page 13.
2. Enter the desired value for the signal level, e.g. 0.32 AU, the delay volume if needed, and select E1 for relay output.
3. Return to the operation screen
4. Move cursor to field Edit Prog, press any numeric key to generate a identifier for you new program. Press ∇ to confirm and to enter the edit mode for the new program.
5. At the time step 0 min the Events-field rests unchanged. At time 5 min program Events 1 0, see Fig. 17 on page 20 and prepare deactivating of the event at 10 min by entering 0 0 again.

Event 1 (relay output) will be enabled from 5 to 10 minutes. Any peaks within the time window exceeding the signal level (0.32 in this example) defined in the SETUP menu leads to generation of a event at the signal outputs.

Peaks exceeding the signal level between 0 and 5 minutes will not produce an event signal.

The value for level chosen in the SETUP menu, will always apply for all programs, it cannot be set individually for single programs.


Fig. 16 Four steps for creating a program

Program Execution

Programs must be activated prior to execution. Enter the desired program number in the field **Load Prog** and confirm by arrow key. The instruments settings are adjusted to all presets of that program. An autozero is performed automatically every time a new program is loaded.

The time table of the program is started by moving the cursor to the **ON** field in the lower left corner of the display. Pressing **1** then starts the program. The display switches to **RUN**, and the elapsed time is displayed in the field Time. Program execution can be interrupted by pressing **0**. It is stopped with the actual values, indicated by a blinking time line. The message **HOLD** is displayed.


Fig. 17 Programming steps for example program no. 7

Execution of the program can either be continued by pressing **1** and returning to **RUN** again or it can be stopped by pressing **0** again. If you stop program execution the instrument will return to the parameters of the start time. These initial parameters are activated also after a program run is finished.

A program can also be started via input START IN on the rear panel of the instrument. In case a trigger signal is applied to this input, the system switches into mode **RUN**.

Creating Program Links

Program links are cross references to existing programs, which can be defined and edited similar to programs themselves. Links cannot reference other links.

In the editing mode links are defined and loaded as the reserved program numbers 21 to 29, where program 21= L1, 22 = L2 etc.

Links have line numbers from **No. 01** up to **No. 99**, every line containing one program, one number for repetitions and a wait-status programmable. The K-2600 UV detector can store up to 9 different links, but the maximum number of programs used in all links is limited to 100. Each program may appear in any frequency in any link.

Programming a link

Enter the number of the desired link in the field Edit Prog, for example 22 for link L2 and press any arrow key to confirm. Part A in Fig. 18 is showing that step, the cursor blinking at **No****. Now press any numeric key to generate a new, continuously numbered line. In this case **No01** will appear.

Then use **▷** to move the cursor to **P**-field and choose the desired program number, e.g. 03. Press again **▷** to put in at the **R**-field the number of runs to be performed with the selected program, e.g. 05.

Finally, in the field **Wait** the values 1 or 0 can be entered. For Wait = 1 the photometer will wait for an external start signal or for an operator pressing 1 on the field **RUN** until this line will be performed. During waiting the photometer is switched to the wait-status. If Wait is set to 0, the program lines of the link are worked out continuously without any pauses.

Part B in the figure is showing that programming step.


Fig. 18 Editing a program link

Running program links

Before a program link can be executed, it has to be loaded. For this purpose enter the desired link number in the Load Prog-field, e.g. 22 for L2 and confirm by pressing any arrow key.

During the execution of a link the LINK menu is available in addition to the common display. In the LINK menu the momentary status of the running link is shown. Move cursor to the **Load Prog** -field and press < to enter this menu page. For more explanations refer to Fig. 11 on page 14.

To return from the LINK menu press >.

Deleting Programs and Program Links

To delete a program or a link, it has to be loaded first. Entering 0 in the Edit Prog and pressing a arrow key will produce the message **Delete this program? Confirm by cursor** in the display. Pressing any arrow key will delete the loaded program or link.

Pressing a numeric key or just waiting for 60 seconds without pressing any key will abort the deletion process.

Wake Up-Program

The program # 30 is reserved for the Wake Up function, enabling the start of the lamp or a wavelength program at a preset time or date. The Wake Up function can be used e.g. for ensuring switch on and preheating of the lamp after switching it off overnight. Editing program # 30 will show the following screen:

If the wavelength range is changed, spectral data in the memory of the K-2600 UV detector are erased.

Scan Autozero

Before scanning spectra the K-2600 UV detector must perform a scan autozero. The scan autozero measures the spectrum of the eluent. These data are stored for correction of the spectra of your analytes with the spectrum of your solvent. Thus, a new scan autozero should be performed each time when the solvent is changed. A scan autozero is performed by pressing the red Auto Zero button for 1 sec. or longer. The autozero is then measured and calculated over the wavelength range set for scanning.

Performing wavelength scans

For acquiring scan data, enter the scan menu by pressing the yellow Scan button once. When pressing the Scan button a second time, a scan is performed immediately.

For taking different scans during a run, please be sure to adjust the scan number between scans. Stored scans are overwritten without any message, when the Scan button is pressed twice.

Output of Scan Data

The stored scan data can be given out via analog outputs and with defined data rates. Time scaling of the data can be adjusted between 1 and 50 steps per second (that is approximately 1-50 nm/s) in the SETUP menu of the photometer, see section "SCAN OUT" on page 12. Low output speed settings lead to enhanced spectral resolution of the data output.

To release the output of scanned data, move the cursor using \triangleright to the **OUT** field in the scan menu and press Δ or ∇ to start the scan output. Make sure that the desired scan number is chosen. Since every scan is kept in internal storage it can be edited or given out several times. Scans are deleted by overwriting with new scan data and in case the wavelength range was changed.

Connecting other Instruments to the K-2600 UV detector

Using the Event and Remote Sockets

On the rear panel side of the K-2600 UV detector two electrical connector sockets are located, see [5.8] and [5.9] in Fig. 3 on page 9, which serve to send or receive signals from other instruments. The event and remote interfaces are controlled by programs or the SETUP menu.

Please avoid touching the electrical contacts of the socket lines. Electrostatic discharges when touching the contacts could damage the electronics of the device.

Assembling Plug Strips

For the electrical connections mentioned in the next chapters plug strips with 2, 3 or 4 connectors, enclosed in the accessories, may be mounted as follows:


Fig. 21 Assembling plug strips

SOP 4 Assembling Plug Strips

1. Insert the rounded end of the lever latch into the square opening of the selected connector of the plug strip
2. Press the catch down as indicated by arrow
3. Insert the uninsulated end of the cable into the opening under the catch
4. Release the catch and remove the lever latch from the plug

The cable is now firmly anchored in the plug strip.

Event Socket

The photometer provides two control outputs. Event 1 (E1) is a relay output, Event 2 (E2) a TTL (Transistor-Transistor-Logic) output. All voltages have to be mounted between **GROUND** and the corresponding event. The conditions 0, 1 or are programmed with the numeric keys 0, 1 or \perp .

The event socket can be used e.g. in combination with the fraction collecting function of the photometer. A KNAUER 16-port switching valve or a fraction collector can be connected to the relay output. The relay will close when the corresponding event E1 is set to 1 in the method and if the signal level chosen for fraction collection in the Setup menu is exceeded, thus producing a short circuit that will switch the valve to the next position.

Fig. 22 Event and remote socket connections

EVENT 1 :	Relay: Possible conditions are :	
	Contacts open	0
	Contacts closed	1
	Pulse: (contacts closed for 500 ms)	
EVENT 2 :	TTL (max. consumption 40 mA, $U_{CEmax} = 24$ V). Possible conditions are:	
	Low (< 500 mV)	0
	High (> 2 V)	1
	Pulse (> 2 V for 500 ms)	
ERROR IN:	when receiving a 0 V error signal, the message Error signal was detected appears on the display, programs are stopped and the lamp is switched off.	
ERROR OUT:	an open collector output remains active so long as an error is detected, e.g. the lamp does not start.	
START IN:	if an electrical signal is received, a program or link is started, either by a low signal or by a short circuit	
START OUT:	an open collector output is active for 500 ms, when a program or link is started.	
EXT λ and AUTOZERO	allow for external analog control of the photometer, if the option ANALOG is chosen in the SETUP menu die Option gewählt wurde, see section „CONTROL“ on page 13.	

Remote Socket

The remote socket on the rear side of the K-2600 UV detector serves for sending and receiving start signals or error signals to or from other instruments. For example an injection valve or an autosampler giving start signals can be connected to the START IN input. All voltages have to be mounted between GROUND and the corresponding output.

RS 232 Serial Interfaces

The two RS 232 serial interfaces on the rear side of the device, see [5.4] in Fig. 3 on page 9, enable digital data transfer between the K-2600 UV detector and compatible instruments such as pumps or a PC, equipped with EuroChrom software. All the devices are communicating with each other and form an integrated network.

The KNAUER Net

The photometer can be combined with further instruments from KNAUER in a digital ring connection scheme, enabling control of these instruments by the KNAUER EuroChrom software. The external software control is subject of the following chapter.

Fig. 23 Network of three instruments

All devices are connected in a token ring in the following way: The first RS 232 interface (the upper one) of the first device is connected to the second RS 232 interface (the lower one) of the next device and so on. The ring is closed by connecting the RS 232 interface of the last device to the free interface of the first one.

List of display messages

The following list shows display messages possibly appearing on your K-2600 UV detector in the order of appearance in this manual.

HPLC SPECTRALP HOTOMETER	7
Start time is fixed	17
Delete this line? Confirm by cursor	17
Delete this program? Confirm by cursor	21
Switch lamp OFF? Confirm by cursor	22
Error signal was detected	25
KNAUER NET	27

Software Control of the K-2600 UV detector

The full capabilities of the K-2600 UV detector are accessed under operation with one of the software packages EuroChrom® for Windows or ChromGate®.

Fig. 24 HPLC-software packages

In combination with this software, the K-2600 UV detector can measure with up to four wavelengths simultaneously and perform an (almost) unlimited number of scans over any desired wavelength range in every run. Trouble free detection and instrument control is achieved by the fully digitally operating KNAUER net.

This chapter will explain the basic knowledge needed for working with the K-2600 UV detector under EuroChrom software control. For more detailed informations concerning as well the features of the software as running it, please consider the users' handbook of the software.

Preparing for Software Control

Installing the KNAUER Net

The instruments of the „Maxi“ series, i.e. the K-2600 UV detector and the pumps K-1000 or K-1001, can be connected to the computer in digital communication network, the KNAUER net, see chapter „RS 232 Serial Interfaces“ on page 25.

The serial port of the computer is connected to one of the serial ports on the first instrument, e.g. the photometer, with an Y-shaped cable delivered with the software.

The second port of the photometer is connected to one serial port of the next device, e.g. the pump. The ring is closed by connecting the remaining port of the pump to the second plug of the Y-shaped cable.


Fig. 25 Network of different instruments in the KNAUER net

Configuring the KNAUER Net

In the SETUP menu of the K-2600 UV detector the option external control: KNAUER NET must be selected, see section „CONTROL“ on page 13.


Fig. 26 Configuring the KNAUER net inside the SETUP menu

Simple Maintenance

Control of the lamp's functionality

The deuterium lamp used with the photometer K-2600 has an extended life time to ensure long-time functionality and reliable measurements with low noise and baseline drift as well as high sensitivity. The actual using time of the lamp depends on different factors, like the number of lamp starts, the average burning time and your requirements concerning noise and sensitivity.

To check the functionality of the lamp, the two intensity values **sig** and **ref** to be found in the menu Signal in the SETUP menu „SIGNAL“ on page 12 provide helpful information. The ref value refers to the intensity of the light measured in the reference channel and can be used for checking the quality of the lamp.

On delivery of the photometer, the sig.-value is set to the range of 0.39 to 0.78 and the ref.-value to 0.2 to 0.4, measured with the dummy cell at a wavelength of 240 nm.

We recommend to check the ref value at regular intervals under the conditions mentioned above (dummy cell, $\lambda = 240$ nm). This applies especially in case higher noise levels or decreased sensitivity are observed on working with the photometer K-2600. In these observations coincide with a **ref** value of approx. 0.1 or less, a new deuterium lamp should be installed, see SOP "Changing the lamp".

Changing the lamp

Remove the power plug before opening. Please let the lamp cool down for at least 15 minutes after switching it off.

SOP 5

Changing the lamp

1. Unscrew the housing and remove it by lifting.
2. The lamp is located in a black housing on the right side approximately in the middle of the instrument.
3. Pull the plug of the old lamp, located to the right in the front of the instrument. Unscrew the two screws in the lamp socket with a screw driver and remove the whole lamp from the instrument including its cables.
4. Insert the new lamp assuring that it is correctly seated in the guiding slot and screw it in securely.

Do not touch the glass of the lamp. Should you touch it accidentally, clean it thoroughly with a lint free cloth and i-Propanol.

5. Insert the plug of the lamp. The deuterium lamp has a 3-pole plug, the halogen lamp a 2-pole one. Such it is impossible to plug a lamp into a wrong socket.
6. If needed, adjust the lamp type in the SETUP menu, window LAMP. In this window it is also possible to reset the working hours of the lamp as described in „LAMP“ on page 11.
7. Put the lid back on from the top of the instrument and fasten the screws of the cover.

Cleaning the flow cell

Noisy baselines and low sensitivities may be due to a dirty flow cell. This may also be indicated by a low value for **sig** in the signal menu when flushing the cell with pure solvent. All flow cells can easily be disassembled for cleaning the lenses.


Analytical flow cells

Fig. 27 Cross-section drawing of an analytical flow cell

SOP 6 Cleaning an analytical flow cell

This instruction applies to the analytical flow cells A4061, A4042, and A4045.

1. Unscrew the outer threads with the 3 mm hexagonal spanner enclosed in the flow cell's delivery.
2. Remove the black gasket that carries the lenses with a pair of tweezers or by gently tapping it on a clean surface. The lens is embedded in the gasket and sealed against the flow path with a PTFE seal. This seal should be changed every time when disassembling the flow cell.
3. Take out the lenses and clean them by wiping them with a soft cloth or with an appropriate solvent in an ultrasonic bath. Be careful not to touch the clean lenses with the fingers.
4. Reassemble the cell in the reverse manner, making sure that the PTFE seal does not block the light path.
5. Tighten the outer threads carefully with the spanner in order not to damage the lenses.


Preparative flow cells

Fig. 28 Cross-section drawing of a preparative flow cell

SOP 7 Cleaning a preparative flow cell

This instruction applies to the preparative flow cells A4066, A4067, A4068 and A4069.

The preparative flow cells have a rod shaped light guide instead of the concave lens in the analytical cells.

1. Unscrew the outer thread with a hexagonal spanner
2. Take out the stainless steel plain washer and the PEEK distance plate (not present in A4069).
3. Get hold of the gasket holder with the light guide by grasping it with a pair of tweezers, using the indentations on the outer side of the gasket holder.

4. Push out the light guide and strip the PTFE sealing ring in order to clean the lens.
5. Reassemble the cell in reverse order. Use a new PTFE sealing ring after every disassembly to ensure the consistence of the flow cell.

Adjusting the path length of the preparative flow cells

SOP 8 Changing the path length

This instruction applies to flow cells A4066, A4067 and A4068 with 1/8" and 1/4" connectors. Path lengths can be adjusted to 2, 1.25 and 0.5 mm. On delivery the path length is set to 2 mm. To reduce the path length to 1.25 or 0.5 mm, follow the instructions given.

1. Unscrew the outer thread with the 3 mm hexagonal spanner.
2. Take out the stainless steel plain washer and the PEEK distance plate.
3. Remove the distance plate, put back the stainless steel washer and refasten the thread carefully.

With the distance plate removed, the rod shaped light guide is pushed further into the flow cell (0.75 mm), thus resulting in a reduced path length of 1.25 mm. To reduce the path length further to 0.5 mm, follow the same procedure on the other side of the cell.

To extend the path length again in steps of 0.75 mm, the PEEK distance plates have to be inserted again.

4. Loosen the outer thread, remove the stainless steel plain washers and take out the gasket holder by using a pair of tweezers.
5. Push the light guide approximately 1 mm to the outside to enlarge the path length. Use a clean cloth and do not touch the light guide with the fingers.
6. Put the gasket holder back into the cell.
7. Insert the PEEK distance plate again and then the plain washer.
8. Fasten the outer thread carefully.

When fastening the threads, the rod shaped light guide is pushed back in the correct position inside the cell. Inserting a distance plate thus enlarges the path length for 0.75 mm. It is not necessary to change the PTFE sealing ring when adjusting the path length.

Flow cells with fiber optical connectors

Changing the path length of flow cells with fiber optical connectors can be done in the same way. In these cells the screwthread is replaced by a special adapter (Fig. 29). They additionally contain a lens to focus the light onto the optical fibers.


Fig. 29 Fiber optical connector of a preparative flow cell

Spare parts and accessories

Standard delivery

WellChrom K-2600 UV-Detector without cell

Operation Manual

Power Supply Cable 230 V

RS232 Cable

2 Integrator Cables

Plug Strips with Connectors

Spare Parts

A4071	Deuterium lamp
A4072	Halogen lamp
A4073	Halogen lamp with socket (needed for first installation of a halogen lamp)
M1642	Power supply cable
A0884	RS-232 connection cable
A0755	Y-shaped RS-232 cable for connection with a PC
A1402	Set of plug strips (3x4; 2x3; 1x2 connections) including lever latch
A1404	Set of signal cables
A1467	Flat line cable (10 polig, 1,5m)
G1023	Integrator cable (1 x CINCH)
M1588	Analog connector cable, 2m (2 x CINCH)

Available Flow cells for your K-2600 UV detector

All flow cells are also available equipped with fiber optical connectors for the use with the fiber optics version of the photometer K-2600.

Analytical Cells

Order No. Cell type	Layer Thickness (mm); Connector	ID channel (mm)	Volume (μ L)	Material	Flow Range (mL/min)	Maximum Pressure (bar)
A4061	10 mm; 1/16"	1,1	10	stainless steel, with heat exchanger	20	300
A4042	3 mm; 1/16"	1,0	2	stainless steel	50	300
A4045	3 mm; 1/16"	1,0	2	PEEK	50	30

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Preparative Cells

N0076	rod shaped lens
P2644	gasket holder
P2645	distance plate, PEEK
P2646	plain washer, stainless steel
P7003	sealing ring, PTFE

Technical Data

Wavelengths	190 - 740 nm, $\Delta\lambda$ 8 nm
Simult. Wavelengths	4 with digital control, 2 with usage of analog outputs („stand alone“ mode)
Lamp	Deuterium (standard), Halogen (optionally)
Edge filter	370 nm
Wavelengths accuracy	± 1 nm
Scan function	on flow
Range of measurement	0-4 AU
Sensitivity	2×10^{-5} AU at 240 nm and time constant 1.0 s
Noise	1 λ 1×10^{-5} AU at 240 nm, 1.0 s; 2 λ 2×10^{-5} AU at 240/300nm, 1.0 s
Baseline drift	1 λ 10×10^{-5} AU/h at 240 nm, 1.0s; 2 λ 15×10^{-5} AU/h at 240/300nm, 1.0 s
Time constants	0.1/ 0.2 / 0.5 / 1.0 / 2.0 / 5.0 / 10.0 s
Scalable integrator out	$\pm 0,1$ V / $\pm 1,0$ V / ± 10 V 16 steps
Autozero Range	Full Scale
Display	2 x 24 characters
Control	2 RS232 interfaces, event outputs (TTL, OC, relays), remote connectors, analog input
Miscellaneous	Storage of up to 20 programs,
Weight / Dimensions	5.5 kg, 160 x 185 x 240 mm (W x H x D)
GLP Support	Detailed trace report with operating hours of total, lamp, servo motor; number of lamp ignitions, comprehensive service information

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Manufacturer's address:

Hegauer Weg 38
14163 Berlin, Deutschland

WellChrom Spectro Photometer K-2600, Order Number **A 4270**
complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG);
EN 61010 – 1 (1993)
- Engineering Guidelines (89/392/EWG)
- EMV Ordinance (89/336/EWG)
EN 50081 – 1 (1992)
EN 55011 (1991) Class B
EN 55022 (1987) Class B
EN 50082 – 1 (1992)
IEC 801 – 2 (1984),
includes IEC 41 B (sec) 81 (1992)
IEC 801 – 3 (1984)
IEC 801 – 4 (1988)

The product was tested in a typical configuration.

Berlin April, 10th 2000

Bernward Rittgerodt (Managing Director)

The CE Shield is attached to the rear of the instrument.

Guarantee statement

The guarantee period of the WellChrom Spectro Photometer K-2501 is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Claims under this guarantee are valid only if the enclosed guarantee card is returned to us at the address shown below within 14 days of receipt of the instrument. Defective detectors should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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		Load Prog	Edit Prog Clear 0	Events: 2 off 0 on 1	Signal 1 [au]	Signal 2 [au]									
A	Link View	<table border="0" style="width: 100%; text-align: center;"> <tr> <td style="width: 15%;">01 03</td> <td style="width: 15%;">--</td> <td style="width: 15%;">---</td> <td style="width: 15%;">---</td> </tr> <tr> <td>ON</td> <td>*** *</td> <td>---</td> <td>---</td> </tr> </table>				01 03	--	---	---	ON	*** *	---	---	Setup	GLP
01 03	--	---	---												
ON	*** *	---	---												
		Hold 0 Run 1	Time [min]	λ 1 [nm]	λ 2 [nm]										
		Load Prog	Edit Prog Clear 0	Events: 2 off 0 on 1	Signal 1 [au]	Signal 2 [au]									
B	Link View	<table border="0" style="width: 100%; text-align: center;"> <tr> <td style="width: 15%;">01 03</td> <td style="width: 15%;">00</td> <td style="width: 15%;">---</td> <td style="width: 15%;">---</td> </tr> <tr> <td>ON</td> <td>000.0</td> <td>220</td> <td>254</td> </tr> </table>				01 03	00	---	---	ON	000.0	220	254	Setup	GLP
01 03	00	---	---												
ON	000.0	220	254												
		Hold 0 Run 1	Time [min]	λ 1 [nm]	λ 2 [nm]										
		Load Prog	Edit Prog Clear 0	Events: 2 off 0 on 1	Signal 1 [au]	Signal 2 [au]									
C	Link View	<table border="0" style="width: 100%; text-align: center;"> <tr> <td style="width: 15%;">01 03</td> <td style="width: 15%;">00</td> <td style="width: 15%;">---</td> <td style="width: 15%;">---</td> </tr> <tr> <td>ON</td> <td>010.0</td> <td>240</td> <td>254</td> </tr> </table>				01 03	00	---	---	ON	010.0	240	254	Setup	GLP
01 03	00	---	---												
ON	010.0	240	254												
		Hold 0 Run 1	Time [min]	λ 1 [nm]	λ 2 [nm]										
		Load Prog	Edit Prog Clear 0	Events: 2 off 0 on 1	Signal 1 [au]	Signal 2 [au]									
D	Link View	<table border="0" style="width: 100%; text-align: center;"> <tr> <td style="width: 15%;">01 03</td> <td style="width: 15%;">00</td> <td style="width: 15%;">---</td> <td style="width: 15%;">---</td> </tr> <tr> <td>ON</td> <td>020.0</td> <td>280</td> <td>254</td> </tr> </table>				01 03	00	---	---	ON	020.0	280	254	Setup	GLP
01 03	00	---	---												
ON	020.0	280	254												
		Hold 0 Run 1	Time [min]	λ 1 [nm]	λ 2 [nm]										

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Dr. Ing. Herbert Knauer GmbH