**WellChrom**  
**Pneumatic-Pump K-1900**  
**Pneumatische -Pumpe K-1900**  
**Manual / Handbuch**

V7027, 03 / 2007

## Contents

General Description	3
Using this Manual	4
Conventions in this Manual	4
SOP's in this manual	4
Technical Data of Pneumatic Pump K-1900	5
Setup of the Pneumatic Pump K-1900	5
Unpacking	5
Front View of the Pneumatic Pump K-1900	6
Connecting the Pneumatic Pump K-1900	7
Compressed Air Connections	7
Pump Heads and Hydraulic Connections	7
Identifying installed Pump Heads	7
Pump Heads and Hydraulic Connections	8
100, 250 and 500 ml Head	8
1000 ml Head	8
Setting up the K-1900	9
Piston Backflushing	9
Air Pressure Influences	10
Calibration for Constant Flow	11
Easy Maintenance	12
Changing and Cleaning the Pump Head	12
Packing List	16
Spare Parts and Accessories	17
Warranty statement	33
Declaration of conformity	34

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## General Description

The Pneumatic Pump K-1900 has been developed for solvent delivery within the following fields:

- Process chemistry
- Column packing
- Preparative chromatography.

It is powered by pressurized air. It can deliver liquids with an adjustable flow rate. Since no electrical connection is required, the Pneumatic Pump K-1900 is suitable for operation in explosion protected environments.

The new Pneumatic Pump K-1900 has been developed for constant pressure applications. It is suitable for packing high performance columns (microbore, analytical and preparative) up to 1000 bar. The flow is highly reproducible and only little maintenance efforts are required.

It features a double-piston design and can be configured with four different pump heads to meet a broad range of requirements in flow and pressure ranges:

- 100 ml pump head (0.1 – 99.9 ml/min), stainless steel, up to 1000 bar
- 250 ml pump head (0.1 – 249.9 ml/min), stainless steel, up to 670 bar
- 500 ml pump head (0.1 – 499.9 ml/min), stainless steel, up to 225 bar
- 1000 ml pump head (0.1 – 999.9 ml/min), stainless steel, up to 160 bar

The names of the pump heads have historical reasons. You can also operate any of the mentioned pump heads with our Preparative HPLC Pump K-1800. In the case of the Pneumatic Pump K-1900 the piston strokes operate with a higher frequency and therefore flow rates can be much higher than specified on the pump head.

The characteristics of flow and liquid pressure for the K-1900 with the different pump heads are listed in Table 2, „Technical data of individual pump heads“ on page 5.

Pump heads can be interchanged by loosening only four screws. Experienced users can disassemble the heads for simple maintenance tasks like cleaning, changing piston seals, etc.

A cooling device (heat exchanger) to be mounted on the pump head is available. It allows to operate the Pneumatic Pump K-1900 with supercritical solvents as CO<sub>2</sub> for Extraction and Supercritical Fluid Chromatography.

The power unit of the pump consists of a pneumatic piston which is set in motion by compressed air. The pneumatic power will be transmitted to the piston rods by the power transmit lever.

## Using this Manual

### Conventions in this Manual

**Special Warnings** are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates **helpful advice's**.

**Important Hints** are marked by the marginal hand symbol.

### SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks in the operation of your WellChrom Pneumatic Pump K-1900. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and be copied, applied signed, and filed in order to document the performance of the instrument.

**Please operate the instrument and all accessories according to instructions and SOP's in this manual. This ensures proper results and longevity of your equipment.**

Table 1 *SOP's in this manual*

SOP 1	Preparing your K-1900 for piston backflushing	9
SOP 2	Start of Solvent Delivery	10
SOP 3	Removing the Pump Head	12
SOP 4	Removing and Checking Piston Rods	12
SOP 5	Disassembling the Pump Head	13
SOP 6	Assembling the Pump Head and Exchanging Piston Seals	14
SOP 7	Installing the Pump Head	14
SOP 8	Cleaning and Replacing Check Valves	15

## Technical Data of Pneumatic Pump K-1900

Table 2 Technical data of individual pump heads

Characteristics	Pump Head			
	100 ml	250 ml	500ml	1000 ml
Maximum Flow Rate (ml/min) at Minimum Back Pressure	300	500	1000	2000
Maximum Back Pressure (bar) for Solvent Delivery	1000	670	250	160
Maximum Allowed Air Pressure (bar)	4.2 *)	7	7	7
Piston Stroke (mm)	9.5	9.5	9,5	9.5
Stroke Volume per Cycle (ml)	0.6	1.35	3,05	5.4
Diameter (mm)	6.35	9.525	14,29	19.05
Pneumatic Transformation Ratio	1:250	1:96	1:47	1:23

\*) See warning below.

**With the 100 ml pump head the external air pressure must under no circumstances never exceed 4.2 bar. Exceeding an external air pressure of 4.2 bar will lead to damage of pump head and pump mechanics. This damage is not covered by guarantee.**

## Setup of the Pneumatic Pump K-1900

### Unpacking

Unpack the instrument and check pump and accessories for any damage due to shipping. Place any claims referring to damage at the transportation company responsible for shipping.


Please use the packing list and make sure that the delivery of the Pneumatic Pump K-1900 is complete. If you have any questions regarding complete delivery, please call our service department as soon as possible.

Please remove the protection foil from pump outlet.

For the complete list of spare parts available please refer to chapter „Spare Parts and Accessories“ on page 17.


## Front View of the Pneumatic Pump K-1900

All operating elements and hydraulic connections are located on the front of the pump, see Fig. 1 and Fig. 2.

- 1 Pump head
- 2 Flushing Valve
- 3 Eluent pressure gauge
- 4 Air pressure gauge
- 5 Regulator for compressed air

Fig. 1 Front View of the Pneumatic Pump K-1900

- |                                      |                                 |
|--------------------------------------|---------------------------------|
| <b>1 Pump Head</b>                   | <b>2 Flushing Valve</b>         |
| 1.1 Piston backflushing capillaries  | 2.1 Inlets to flushing valve    |
| 1.2 Eluent inlets                    | 2.2 De-aeration capillary 1/16" |
| 1.3 Head set screws                  | 2.3 De-aeration screw           |
| 1.4 Label with type of pump head     | 2.4 Eluent outlet to column     |
| 1.5 Inlets valve housing             |                                 |
| 1.6 Outlets valve housing            |                                 |
| 1.7 Eluent outlets to flushing valve |                                 |
| 1.8 Connection capillaries           |                                 |

Fig. 2 Pump head and flushing valve

## Connecting the Pneumatic Pump K-1900

### Compressed Air Connections

Make sure that your compressed air system to be used with the Pneumatic Pump K-1900 is able to maintain pressures up to 7 bar. Without back pressure the air consumption would be 80 L/min. If the Pneumatic Pump K-1900 works against pressure the consumption of air is much lower depending on the liquid pressure.

The connector to the compressed air is located on the lower rear panel of the pump. Connect the coupler for the pressure tube to the coupler socket. Push the pressure tube (not delivered) from the compressed air supply over the coupler and secure it with an appropriate clamp. Connect the coupler socket to the compressed air connector of the pump (see Fig. 2).

To disconnect from the compressed air supply: Hold the brass part of the coupler with one hand and pull the blue part of the coupler socket away from the pump with the other hand.


Fig. 3 Air pressure connection of the K-1900

### Pump Heads and Hydraulic Connections

**The pump should be rinsed with a methanol-water mixture (1:9 v/v) before first usage. Make sure that all hydraulic connections are suitable for the system pressure and flow rate of your liquid transportation system.**

#### Identifying installed Pump Heads

Maximum flow rate, 100, 250, 500, or 1000 ml/min, is indicated on a label on the front of the pump head. These values are valid for the HPLC-Pump K-1800. Working with the pneumatic pump K-1900 the maximum flow rates are approximately two times higher.

Fig. 4 Pump head label (1000ml head)

## Pump Heads and Hydraulic Connections

All operating elements and hydraulic connections are located on the front of the pump, see Fig. 1 and Fig. 2.

**Make sure that all hydraulic connections are suited for the system pressure and flow rate of your HPLC System. During initial set-up procedures, and also later when changing solvent, pay strict attention to the miscibility of solvents to be used one after another, otherwise the function of the check valves may be adversely affected (see Fig. 14 on page 15). By substitutions non-miscible solvents such as methanol and hexane, rinse after using methanol with an intermediate solvent that is miscible with both solvents in all proportion, e.g. isopropanol.**

### 100, 250 and 500 ml Head

Solvent tubings are connected to the pump head according to Fig. 5. Make sure that the flat side of the ferrule is directed to the pump head and that the fitting screw is fastened hand tight. According to the local conditions the tubings can be connected either to the inlets valve housing (A) or to the eluent inlet (B).


Fig. 5 100, 250 and 500 ml pump head: hydraulic connections

### 1000 ml Head

Push the tubes directly over the connectors, and secure them by means of appropriate clamps.

Fig. 6 1000 ml pump head: hydraulic connections

## Setting up the K-1900

**Never run the Pneumatic Pump K-1900 without liquid in the pump head or in the piston backflushing compartment. For piston backflushing, follow the instructions below. Operation with-out liquid may lead to damage of the pump seals in short time.**

### Piston Backflushing

Backflushing the piston removes traces of salt and other decontaminants from the backside of the pistons.

Especially if you use saline solvents or buffer solution we highly recommend that you continuously backflush in order to prevent crystallization which can damage your piston seals.

For best protection and error free long time operation of the pump head it is recommended always to backflush the piston with an appropriate solvent (mostly water) according to the following SOP.

### SOP 1 Preparing your K-1900 for piston backflushing

This SOP applies to the Pneumatic Pump K-1900.

1. Push a 1/16" ID tube over both piston backflushing capillaries.


Fig. 7 *Piston backflushing capillaries*

2. Place one low end of the tubing in a flask.
3. Fill the priming syringe with the backflushing liquid.
4. Connect the syringe with the other tubing.
5. Press liquid through the pump head, until it flows without any air bubbles into the container.
6. Remove tubings from capillary opening.

If you want operation with continuous backflushing you can attach two containers of rinsing liquid instead of the priming syringe. The containers should be positioned that one container is located higher than the other, thus ensuring liquid flow through the pump head without any assistance.

When not backflushing continuously, we recommend to connect both capillary openings with a tubing to prevent vaporization of solvents and drying out of the piston chamber.

## SOP 2 Start of Solvent Delivery

This SOP applies to the Pneumatic Pump K-1900.

1. Screw a blind fitting to the "Eluent outlet to column", Pos. {2.4} in Fig. 2 and open the "De-aeration screw", Pos. {2.3} for opening flushing valve.
2. Connect the Luer-Lock syringe with the "De-aeration capillary", Pos. {2.2} in Fig. 2 with the help of a short PTFE tube (ID 1,5 mm/ OD 2,1 mm). The pump head needs to be filled with liquid by pulling eluent into the Luer-Lock syringe until liquid enters the syringe.
3. Close the "De-aeration screw" and exchange the blind fitting with a regular 1/8" fittings to connect this outlet with an injection system, a column packing device, a chemical reactor or something else.
4. Pull the Regulator for compressed air, see Fig. 2 for loosening up. The pneumatic can be set in operation by turning the Regulator for compressed air to the right (required air pressure increases). The eluent will be transported to flushing valve.
5. The pump should be stopped by turning the Regulator for compressed air to the left (required air pressure decreases towards 0 bar) and any unit (e.g. the column filling device) is connected to the solvent outlet.

### Air Pressure Influences

The resulting eluent pressure increases or decreases as a function of compressed air inlet pressure. Eluent pressure will nearly be constant, if the compressed air inlet pressure remains constant. The *Regulator for compressed air* can be locked by pressing it to wards the instrument.

The delivery volume depends on the required air pressure and the loss of pressure along the column or any device that generates back pressure. You can reach the maximum delivery volume at minimum eluent pressure below 1 bar and the maximum eluent pressure at the minimum flow rate. Fig. 8 up to Fig. 10 show the delivery volume as a function of the eluent pressure for the 100, 250 and 1000 ml pump head, respectively. In those figures you can see the limit values for the liquid pressure and the flow rate corresponding to Table 2, "Technical data of individual pump heads," on page 5 and the maximum allowed air pressure for each pump head. By the constant required air pressure the flow rate increases if the liquid pressure decreases. The Figures should help to choose the pump head for your application.


Fig. 8 Interdependence of delivery volume and eluent pressure for the 100 ml pump head


Fig. 9 Interdependence of delivery volume and eluent pressure for the 250 ml pump head


Fig. 10 Interdependence of delivery volume and eluent pressure for the 1000 ml pump head

## Calibration for Constant Flow

If you want to use the pneumatic pump as a delivery pump with constant delivery volume, you need to calibrate it.

Fig. 11 shows typical calibration curves for two different columns with different pressure losses along the column for the 250 ml pump head. In this figure column 1 has a lower back pressure than column 2 at same flow rate, so it is possible to obtain higher flow rates with the same required air pressure on the column 1 compared to the column 2.

For the calibration you have to measure the flow rate at different required air pressures from 1 up to 7 bar (up to 4.2 bar with the 100 ml pump head)


Fig. 11 Typical calibration curves of the Pneumatic Pump K-1900 for two different columns (250 ml pump head)

## Easy Maintenance

### Changing and Cleaning the Pump Head

#### SOP 3 Removing the Pump Head

This SOP refers to the Pneumatic Pump K–1900.

For exchange of the pump head, or for disassembling it in order to clean valves or replace piston seals, piston rods etc., the pump head can easily be removed.

1. Clean the pump head with a suitable cleaning reagent and then with distilled water according to chapter „Start of Solvent Delivery“ on page 10, step 1 and 2.

**If organic solvents remain in the pump head, danger of skin irritation may exist.**

2. Remove the two solvent tubings from the „Eluent inlets“, Pos. {1.2} see Fig. 2 on page 6.
3. Loosen the two „Eluent outlets to transducer“, Pos. {1.7} and the two „Inlets to transducer“, Pos. {2.1}.
4. Remove the two „Connection capillaries“, Pos. {1.8}, beginning with the upper „n“-shaped one.
5. Loosen the „Inlets valve housing“, Pos. {1.5} as well as the fittings at „Outlets valve housing“, Pos. {1.6}.
6. Using a hexagonal spanner no. 4 (4 mm), loosen just two diagonally opposed „Head set screws“, Pos. {1.3}. Remove the screws.
7. Carefully loosen the two remaining screws, alternating from one to the other, approx. half a turn. This prevents the pump head from tilting and becoming damaged. Once the spring tension has been reduced, hold the pump head firmly in one hand while removing the screws completely with the other hand.
8. Carefully remove the pump head.

#### SOP 4 Removing and Checking Piston Rods

This SOP refers to the following pump heads:

100 ml stainless steel, KNAUER order number A 4029  
 250 ml stainless steel, KNAUER order number A 4021  
 500 ml stainless steel, KNAUER order number A 4038  
 1000 ml stainless steel, KNAUER order number A 4022

1. Remove Pump head as described in SOP „Removing the Pump Head“.

**If you only wish to check the piston rods, you do not need to disassemble the pump head any further.**

2. The „Piston rods“, Pos. {3.17} in Fig. 12 on page 13 may be removed using pliers. Take the tip of the piston using the pliers, and pull it out carefully in a straight line.

**If the rods are broken, you must check the entire pump head for damage.**

## SOP 5 Disassembling the Pump Head

This SOP refers to the pump heads 100 ml, 250 ml, 500 ml, and 1000 ml.


Fig. 12 Explosion diagram of pump heads

1. Apply SOP 4 "Removing and Checking Piston Rods".
2. Loosen the two „Retaining plate screws“, Pos. {3.19} half a turn, alternating from one to the other to avoid damaging the „Retaining plate“.
3. Because the two screws are very tight, it may be helpful to either clamp the pump head or to press one of its side surfaces against a table with one hand while loosening the screws.
4. Unscrew the two screws strictly alternating due to strong force of the „Springs“, Pos. {3.15} behind the plate, and remove them.
5. Remove the „Retaining plate“, Pos. {3.18}.
6. You can now remove the „Guides for springs“, Pos. {3.16}, the „Springs“, Pos. {3.15} as well as the two „Distance plates“, Pos. {3.13}.
7. Using the SW 10 spanner, loosen and remove the „Spacing bolts“, Pos. {3.14}. These bolts are seated very tightly. Follow the tip given in step 3.
8. Lift and remove the „Pump head housing, basic plate“, Pos. {3.B}.
9. To remove the piston seals, Pos. {3.11} in the basic plate {3.B} and in the front block {3.A}, use an appropriate tool, e.g. a screw driver and carefully lever the four piston seals.

## SOP 6 Assembling the Pump Head and Exchanging Piston Seals

This SOP refers to the pump heads 100 ml, 250 ml, 500 ml, and 1000 ml.

All positions of components refer to Fig. 12 „Explosion Diagram and list of parts of 100 ml, 250 ml, and 1000 ml pump heads“ on page 13.

1. It is recommended to exchange all four piston seals located in the pump heads after disassembling the pump head.

Fig. 13 Open side of piston seal

2. With the open side facing downwards, carefully press the „Piston seals“, Pos. {3.11} into the „Pump head housing, front block“, Pos. {3.A} and „Pump head housing, basic plate“, Pos. {3.B}, making sure to keep them straight.
3. Insert the O-ring sealings part {3.12}.
4. Put the two housing parts together, positioning the O-ring side of part {3.B} against part {3.A}. Be sure that the backflushing capillaries are positioned on the top side of the head - see label on front of part {3.A} for orientation.
5. Insert the two „Spacing bolts“, Pos. {3.14}, screw them in with the hand and tighten them as firmly as before using the SW 10 spanner.
6. Put in the two „Distance plates“, Pos. {3.13}, position the „Springs“, Pos. {3.15} and the „Guides for springs“, Pos. {3.16}.
7. Put the „Retaining plate“, Pos. {3.18} over the „Guides for springs“, Pos. {3.16}, with the conical bore of the plate directed outwards.
8. Insert the two „Retaining plate screws“, Pos. {3.19}, press the plate down with one hand and tighten the screws with a screw driver, strictly alternating due to strong force of the „Springs“, Pos. {3.15}.
9. Carefully insert the piston rods, see „Piston rods“, Pos. {3.17} without bending or quenching the rods.
10. The „Retaining plate screws“, Pos. {3.19} must be tightened that they are seated as securely as before.

## SOP 7 Installing the Pump Head

This SOP refers to the following pump heads:

100 ml stainless steel, KNAUER order number A 4029  
 250 ml stainless steel, KNAUER order number A 4021  
 500 ml stainless steel, KNAUER order number A 4038  
 1000 ml stainless steel, KNAUER order number A 4022

1. Make sure that the pump head is properly assembled.
2. Position the head in a straight line onto the pump housing.
3. Tighten all four „Pump head set screws“, Pos. {3.3} a few turns by hand.
4. Alternating from one to the next, tighten two diagonally opposed screws half a turn at a time, until the pump head is correctly seated.
5. Tighten the two other screws.
6. Make sure that all four „Pump head set screws“, Pos. {3.3} are securely tightened.

7. Install the two „Connection capillaries“, Pos. {3.8}, beginning with the lower „s“-shaped one, and firmly tighten all fittings.

## SOP 8 Cleaning and Replacing Check Valves

This SOP refers to the following pump heads:

- 100 ml stainless steel, KNAUER order number A 4029
- 250 ml stainless steel, KNAUER order number A 4021
- 500 ml stainless steel, KNAUER order number A 4038
- 1000 ml stainless steel, KNAUER order number A 4022

If the check valves become dirty they will no longer open and close correctly. You can remove the check valves for cleaning from the pump head after loosening and removing the „Inlet valve housing“, Pos. {3.5} and „Outlet valve housing“, Pos. {3.6}, and then disassemble them according to the following instructions.

1. Place the valves in a suitable cleaning solution. Use a ultrasonic bath to clean the valve. If malfunction persists, use steps 2...5 of this SOP.
2. Using a knife or similar, carefully remove the valve seal from the housing, see Fig. 23 on page 41.
3. Remove the valve ball out of the housing.
4. Clean the individual parts. We recommend a ultrasonic bath.
5. Assemble the check valves in reverse order, use Figure Fig. 14.


Fig. 14 Components of the check valves of all pump heads available

## Packing List

The delivery of the Pneumatic HPLC Pump K-1900 with the 100 ml, 250 ml, 500 ml, or 1000 ml pump head (order number A 40281, A 40282, A 40285, A 40289) comprising the items according to Table 3.

Table 3 The delivery consists of:

Peaces	Description
<b>Pump basic unit:</b>	
1	User manual
1	Coupler, to attach air pressure tube
1	Coupler socket
1	Blind fitting
<b>Set of tools including:</b>	
1	Double open-end wrench, 13 and 17 mm
1	Double open-end wrench, 8 and 10 mm
1	Double open-end wrench, 1/4" and 5/16"
1	Screw key, 4 mm
1	Recessed head screw No. 2
<b>Syringe set including</b>	
1	Syringe
1	Luer-Lock needle
1	0,1 m PTFE tube (ID 1,5 mm/ OD 2,1 mm)
<b>Eluent outlet set including</b>	
1	Capillary AD 1/8"
2	Fitting screws
2	Ferrules
<b>Set of fittings and adapters (1/8" – 1/16") including</b>	
1	1/8" – 1/16" transforming coupler
1	Fitting screw for capillary AD 1/8" (M8x1)
1	Ferrule for capillary AD 1/8"
1	Fitting screw for capillary AD 1/16"
1	Ferrule for capillary AD 1/16"
<b>Pump head 100 ml, order number A 4029 *)</b>	
1	Pump head 100 ml
4	Tube fitting screws
6	Tube ferrules
2	1 m PTFE tube (ID 3 mm/ OD 4 mm)
1	0,5 m PTFE tube (ID 1,5 mm/ OD 2,1 mm)
<b>Pump head 250 ml, order number A 4021 *)</b>	
1	Pump head 250 ml
4	Tube fitting screws
6	Tube ferrules
2	1 m PTFE tube (ID 3 mm/ OD 4 mm)
1	0,5 m PTFE tube (ID 1,5 mm/ OD 2,1 mm)
<b>Pump head 500 ml, order number A 4038 *)</b>	
1	Pump head 500 ml
4	Tube fitting screws
6	Tube ferrules
2	1 m PTFE tube (ID 3 mm/ OD 4 mm)
1	0,5 m PTFE tube (ID 1,5 mm/ OD 2,1 mm)
<b>Pump head 1000 ml, order number A 4022 *)</b>	
1	Pump head 1000 ml
2	1 m PTFE tube (ID 7 mm/ OD 9 mm)
1	0,5 m PTFE tube (ID 1,5 mm/ OD 2,1 mm)

\*) Please select the appropriate pump head for your special requirements.

## Spare Parts and Accessories

Table 4

	KNAUER Order Number
User Manual	V 7027
<b>Optional Pump Heads</b>	
Pump Head 100 ml incl. accessory kit	A 4029
Pump Head 250 ml incl. accessory kit	A 4021
Pump Head 500 ml incl. accessory kit	A 4025
Pump Head 1000 ml incl. accessory kit	A 4022
<b>Spare Parts for Pump Head 100 ml</b>	
see Explosion Diagram of Pump Heads on page 13	
Check valve {3.10}	A 1122
Piston rod 1/4'' {3.17}	A 0747
Set of Piston seals including 2 piston seals for 1/4'' pistons{3.11} 1 O-shaped sealing {3.12}	A 0746
Set of Hydraulic connections including 10 tube fitting screws and 10 tube ferrules	A 2003
Inlet junction including 2 Inlet valve housings 1 T peace, fitting screw with ferrule 1 left-hand and 1 right-hand connection tube	A 1121
<b>Spare Parts for Pump Head 250 ml</b>	
Check valve {3.10}	A 1122
Piston rod 3/8'' {3.17}	A 1017
Set of Piston seals including 2 piston seals for 3/8'' pistons{3.11} 1 O-shaped sealing {3.12}	A 1046
Set of Hydraulic connections including 10 tube fitting screws and 10 tube ferrules	A 2003
Inlet junction including 2 Inlet valve housings 1 T peace, fitting screw with ferrule 1 left-hand and 1 right-hand connection tube	A 1117
<b>Spare Parts for Pump Head 500 ml</b>	
Check valve {3.10}	A 1080
Piston rod 9/16'' {3.17}	A 1478
Set of Piston seals including 2 piston seals for 9/16'' pistons{3.11} 1 O-shaped sealing {3.12}	A 1479
<b>Spare Parts for Pump Head 1000 ml</b>	
Check valve {3.10}	A 1080
Piston rod 3/4'' {3.17}	A 1015
Set of Piston seals including 2 piston seals for 3/4'' pistons{3.11} 1 O-shaped sealing {3.12}	A 1077

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## Warranty statement

The warranty period of the WellChrom Pneumatic Pump K-1900 is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Defective pumps should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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## Declaration of conformity

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**Manufacturer's name and address:  
Herstellername und -adresse**

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complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)  
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)  
EN 6100-3-2 (10/2006)  
EN 61326-1 (10/2006)

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Dr. Ing. Herbert KNAUER GmbH