# Smartline Sample Preparation Unit 6300

V5003 07/2007

*Smartline*

Wissenschaftliche Gerätebau  
Dr. Ing. Herbert Knauer GmbH  
Hegauer Weg 38  
D - 14163 Berlin, Germany  
Tel.: +49 (0)30 809 727 0  
Fax.: +49 (0)30 801 50 10  
E-Mail: [info@knauer.net](mailto:info@knauer.net)  
Internet: [www.knauer.net](http://www.knauer.net)

---

## Content

Using this Manual .....	2
Conventions.....	2
SOPs in this manual.....	2
General Description of the Sample Preparation Unit 6300 .....	3
Smartline Pump 1000 .....	3
Smartline Assistant 6010.....	4
Unpacking.....	4
Power Supply .....	4
Front View of the Smartline Assistant 6010 .....	4
Rear View of the Smartline Assistant 6010.....	5
Gradient Mixing Block .....	5
Valves and Valve Drives.....	6
Operating the Smartline Assistant 6010.....	7
Electrical Connections.....	7
Eluent Connection with the Smartline Pump 1000.....	8
Capillary Connections .....	9
Coupling of the Sample Preparation Unit 6300 to a LC System .....	11
Function of the Sample Preparation Unit 6300 .....	11
Example Program for the Sample Preparation Unit 6300 .....	12
Delivery List and Accessories.....	13
Technical Data.....	14
Warranty statement .....	28
Declaration of conformity .....	29

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## Using this Manual

### Conventions

**Special Warnings** are indicated by the marginal warning sign and printed in bold letters.

**Important hints** are marked by the marginal hand symbol.

The marginal **light bulb** symbol indicates **helpful advice**.

### SOPs in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks for operating Injection and Switching Valves. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and can be copied, applied, signed, and filed in order to document the performance of the instrument.

**It is very important to follow all instructions and SOPs in this manual when operating the valves. This ensures proper results and longevity of your equipment.**

*Table 1 List of SOPs in this manual*

SOP 1	Mounting Plug Strips	8
SOP 2	Capillary connections using the DYNASEAL system	10
SOP 3	Capillary connections using the UNF system	10
SOP 4	Capillary connections at the Smartline Assistant 6010	11

## General Description of the Sample Preparation Unit 6300

The sample preparation unit 6300 includes the Smartline Pump 1000 and the Smartline Assistant 6010.

Fig. 1 Front view of the sample preparation unit SAMPLE PREPARATION UNIT 6300

It offers the possibility to automate all steps of the conventional solid phase extraction with direct elution onto the HPLC column. As a result, the reliability of analytical data will be increased. While the HPLC part performs the chromatographic separation, the sample preparation module SAMPLE PREPARATION UNIT 6300 prepares the next sample.

### Smartline Pump 1000

The Smartline Pump 1000 is a representative of the youngest generation of HPLC pumps. The pump offers very low residual pulsation, high flow rate precision and accuracy. Its optimized mechanics and piston back flushing pump provide for reliability and a long lifetime. The pump head is available with stainless steel, titanium or ceramic inlays. Thus its various applications also as biocompatible option is possible. The pump can be controlled by the foil key touchpad. It offers the possibility to store up to 20 programs (flow rate, gradient and valve switching). Additionally the pump can be controlled by the ChromGate software as well as over an analogue signal.

**Use and operation of the Smartline Pump 1000 is described in the separate pump manual.**

## Smartline Assistant 6010

The Smartline Assistant 6010 includes two 6-port/3-channel-valves and one gradient mixing block, which completes the pump to a quaternary LPG system with highest gradient- and flow precision. Both valves as injection- and bypass valves are controlled by the pump.

### Unpacking

All our instruments are always packed carefully and safely for transportation. After unpacking, please check the device and accessories thoroughly for any damage that may have occurred during transport. If necessary, put forward any claim for damages to the carrier.

Use the delivery list at page 13 and check that the delivery is complete. Please contact your local dealer or our sales department if you are missing something or if you need support.

Remove the transport protection from the display.

### Power Supply

The Smartline Assistant 6010 is equipped with a universal power supply for input voltages ranging from 90 to 260 Volts AC. Manual switching the supply voltage is not required.

**Make sure that the main power supply is properly grounded and that an appropriate three-pronged power cable is used. Connecting the pump to a faulty power socket can cause damage to the instrument.**

Connect the power cable to the socket on the rear of the instrument (see Fig. 4, Pos. 3), and switch on the device with the ON/OFF key (Pos. 2).

### Front View of the Smartline Assistant 6010

Fig. 2 Front View of the Smartline Assistant 6010, door closed


Fig. 3 Front View of the Smartline Assistant 6010, door opened

- 1 Valve displays
- 2 Valve foil pads
- 3 Gradient mixing block, inputs
- 4 Gradient mixing block, outputs to pump head
- 5 Injection valve
- 6 Bypass valve

### Rear View of the Smartline Assistant 6010


Fig. 4 Rear View of the Smartline Assistant 6010

- 1 Socket for the pump connection cable
- 2 Power switch
- 3 Power supply socket
- 4 RS 232 connectors for the valves
- 5 WAGO remote control connectors for the valves

### Gradient Mixing Block

The gradient mixing block complements the Smartline Pump 1000 that creates a completely quaternary low pressure gradient system. The flow range reaches up to 10 ml/min. A valve block made out of PEEK is controlled via the Smartline Pump 1000 and delivers gradient mixtures with the highest precision over the entire flow range.

The LPG module in the Smartline Assistant 6010 can only be operated in combination with a Smartline Pump 1000. The pump controls all valves of the

gradient mixing block and provides the required electrical power supply. This is realized with a connection cable. The connection cable is included in the accessories for the Smartline Assistant 6010.

The control of the gradient mixing block can be obtained in the stand-alone mode via the Smartline Pump 1000 or with ChromGate software.

## Valves and Valve Drives

The valve drives can be operated manually using the foil pad, (Pos. 2 in Fig. 3). Control is also possible in a stand alone mode via the Smartline Pump 1000 and the remote control strips (Pos. 5 in Fig. 4) or by help of the ChromGate® HPLC software via RS 232 connectors (Pos. 4).

**Valves and valve drives are described in detail in the corresponding manuals (V 7447 valve drives and V 7452 valves).**

Table 1 Allocation of the WAGO remote control strips

<b>GROUND:</b>	Ground connector
<b>HOME:</b>	Sets the valve drive to position 1
<b>FORWARD / LOAD:</b>	Sets the valve drive to the LOAD position
<b>BACKWARD / INJECT:</b>	Sets the valve drive to the INJECT position
<b>OUT:</b>	Control output, low-resistance in position INJECT and high resistance in the position LOAD

All KNAUER Injection and Switching Valves were developed for HPLC use. The valves are available in stainless steel and PEEK. If the PEEK version is used, the sample and eluent only come into contact with the polymer. Thus, the prerequisite for working under inert conditions is met.

In general, the valves are adjusted for pressure up to 300 bar. Upon request, they can be adjusted for pressure up to 400 bar. It is also possible for our service department to readjust the valves.

The Multi-Channel Injection and Switching Valves are equipped with a reed contact and magnet. By switching the valve to the inject position a contact is made. Switching the valve to the load position will reopen the same contact.

**Do not use the valves in dry conditions. Otherwise, it can result in a decrease in the lifetime of the rotor sealing. This is especially true of the motor driven valves.**

**If aqueous salt or buffer solutions are used as eluents, it is strongly recommended to rinse the valve with distilled water and/or methanol after use. This will prevent the formation of crystalline deposits inside the valve.**

Fig. 5 shows the internally connected ports according to the valve position in LOAD or INJECT.


Fig. 5 Internally connected ports according to the valve position

## Operating the Smartline Assistant 6010

The Smartline Assistant 6010 and the Smartline Pump 1000 must be arranged in a tower on top the other.

### Electrical Connections

In order to correctly install the electrical connection, carry out the following directions in the order given below:

1. Plug the included connection cable into the "HPLC Manager" outlet on the rear panel of pump, Fig. 6, Pos. 1.
2. Plug the cable into the "pump" outlet on the rear panel of the Smartline Assistant 6010, Pos. 2.

When the Smartline Pump 1000 is turned on, the yellow LEDs A to D will light up for approximately three seconds and then they will shut off.


Fig. 6 Electrical connections of the Sample Preparation Unit 6300

3. Connect the 14 pin WAGO plug to the remote control strip (Pos. 3) of the pump.
4. Connect the 5 pin WAGO plugs to the remote control strips (Pos. 4 and 5) of the Smartline Assistant 6010.

Generally it is your free choice which of the 5 pin plugs is connected to which socket 4 or 5. Only the allocation of the programmed events will be exchanged, compare Fig. 7 and SOP 9 „Setting the Control Signals“ (Events) on page 19 of the pump manual.


Fig. 7 Relation between event programming and the event terminal strip

### Configuring Plug Strips

If necessary the special valve control cable can be configured as described below.

#### SOP 1 Mounting Plug Strips

Use this SOP for connecting wires to the plug strips.

- 
1. Insert the rounded end of the lever latch into the square opening of the selected connector of the plug strip.
  2. Press the latch down as indicated by the arrow.
  3. Insert the uninsulated end of the cable into the opening under the latch.
  4. Release the latch (1) and remove the lever latch from the plug (2).  
The cable is now firmly anchored in the plug strip.

Fig. 8 Mounting the WAGO plug strips

### Eluent Connection with the Smartline Pump 1000

The required tube connections for the four LPG connections from the LPG block to the Smartline Pump 1000 are cut (0.27, 0.28, 0.29 and 0.30 m) to fit and included in the accessories. The black PTFE bushings and seals must be attached to all of the tubes.

Connect the LPG block "OUT" to the low-pressure mixing block on the pump head using the 0.27 to 0.30 m tubes (see Fig. 9). Close the inputs that are not used with blind fittings.


Fig. 9

LPG tubing connections

**Please make sure that the solvent bushings and the blind fittings are tight. If they are not, the specifications of the LPG cannot be guaranteed.**

Set the desired low-pressure gradient in the pump's setup menu according to SOP 3 Selecting a Gradient Mode and SOP 8 Setting the Solvent Composition of the pump's manual.

## Capillary Connections

In general, the capillary connections to the Injection and Switching Valves can be made with stainless steel or PEEK capillaries. In both cases the capillaries can be fixed with the DYNASEAL connection system. However, stainless steel capillaries provide a longer secure connection with the enclosed UNF bushings.

**To avoid damage do not use steel connections for PEEK valves. Only use polymer capillaries, split-grooved clamping rings, polymer sealing rings, and DYNASEAL bushings. The DYNASEAL bushings should be tightened by hand only.**

All capillary connections must be made with minimal dead volume. Therefore, use the shortest capillaries possible with a small inner diameter.

Fig. 10 shows a schematic diagram of the UNF and DYNASEAL connection system.


Fig. 10 Principle of the UNF and DYNASEAL connection systems

Both systems are available with short and long bushings. The shape of the bushings is not relevant to their function. Only make your selection according to the local geometrical conditions.

**The accessibility of the bushings can often be optimized by alternately inserting the bushing types.**

## SOP 2 Capillary connections using the DYNASEAL system

1. Make sure that the capillary is squarely cut off. If necessary, cut the capillary using the tube cutter (Order No. 0569).
2. Push the bushing, the split-grooved clamping ring, and the sealing ring onto the capillary. Pay attention to the sequence and alignment of the fittings (see Fig. 10).
3. Push the capillary into the valve port as far as possible.
4. Hand tighten the bushing while pushing the capillary towards the port (see Fig. 11).

## SOP 3 Capillary connections using the UNF system

1. Make sure the capillary is squarely cut off. If necessary, cut the capillary using the cutting pliers and deburr it with the capillary cleaning set.
2. Push the bushing and the ferrule onto the capillary. Pay attention to the sequence and alignment of the fittings (see Fig. 10).
3. Push the capillary into the valve port as far as possible.
4. Hand tighten the bushing while pushing the capillary towards the port (see Fig. 10).
5. Slightly tighten the bushing with a hexagonal spanner.


Fig. 11 Seating of capillary connectors

#### SOP 4 Capillary connections at the Smartline Assistant 6010

1. Prepare a capillary connection from port 1 of the injection valve to port 5 of the bypass valve.
2. Prepare a capillary connection from port 4 of the injection valve to the cartridge and from there to port 6 of the bypass valve.
3. Prepare a capillary connection from the pump head outlet (see Fig. 9 on page 9) to the autosampler and from there to port 2 of the injection valve.
4. Connect waste tubings as well to port 3 of the injection valve as to port 2 and port 4 of the bypass valve.
5. Connect your LC pump to port 5 of the injection valve.
6. Prepare a capillary connection from port 6 of the injection valve to the column inlet of your HPLC system.

#### Coupling of the Sample Preparation Unit 6300 to a LC System

LC pump, HPLC column, outlet of the gradient mixing block and one side of the cartridge are connected to the injection valve. The other side of the cartridge is connected to the bypass valve. This offers the possibility to fill all capillaries with new eluent to avoid interferences mixing different eluents. Additionally the cartridge can be dried by a compressor.

The sample can be presented onto the cartridge by means with an autosampler. The autosampler is connected to the gradient mixing block and the injection valve. If higher amounts of sample should be enriched onto the cartridge as it is possible with an autosampler, the sample can be presented directly by one channel of the gradient mixing block.

#### Function of the Sample Preparation Unit 6300

The sample is presented onto the cartridge in LOAD position of both valves (Fig. 12). The used eluents for conditioning, elution of the sample and further cleaning of the cartridge are delivered with the gradient mixing block.


Fig. 12 Eluent flow through the valves in LOAD – LOAD positions

In INJECT position of the injection valve (Fig. 13) the transfer of the enriched analyte to the column will be presented. After some minutes this valve can be switched to LOAD position to start the next sample preparation can be filled with new eluent.


Fig. 13 Eluent flow through the valves in INJECT – LOAD positions

In INJECT position of the bypass valve (Fig. 14) the capillaries can be filled with new eluent. Additionally the cartridge can be tried by a compressor


Fig. 14 Eluent flow through the valves in LOAD – INJECT positions

**A switching of both valves into the INJECT positions has to be avoided. In this situation the LC pump and the compressor would operate against each other.**

### Example Program for the Sample Preparation Unit 6300

The given example includes in the first step the sample application by an autosampler followed by a cartridge washing. By switching the valve to INJECT position the enriched sample will be transferred to the LC system.

To clean the cartridge the injection valve is switched back to LOAD position. The cartridge can be dried using a compressor connected to the bypass being switched in INJECT position. At last for conditioning the cartridge the bypass valve is switched again to LOAD.

The programming is performed on the Smartline Pump 1000.

Time Min	Flow rate ml/min	Channel A Eluent 1	Channel B Eluent 2	Events	Operating Steps
0.00	0.40	100 %		00000000	Sample Application
0.02	0.40	100 %		00000000	Wash Step, Eluent 1
4.02	0.40		100 %	00000010	Analyte Transfer
8.02	0.40		100 %	00000000	Wash Step, Eluent 2
12.00	0.40		100 %	00000001	Start Bypass
14.00	0.40		100 %	00000000	End Bypass
14.02	0.40	100 %		00000000	Conditioning with Eluent 1
18.00	0.40	100 %		00000000	End Program

## Delivery List and Accessories

### The delivery consists of:

- Smartline Pump 1000 with installed pump head
- Smartline Assistant 6010
- Bottle tray
- 2 User manuals

**Smartline Pump 1000 accessories see pump manual.**

### Smartline Assistant 6010 accessories:

- Power supply cable
- RS232 cable
- Connecting cable
- Control cable with 3 WAGO strips
- Eluent inlet set
  - 1 fitting screw 1/8"
  - 2 ferrules 1/8"
  - 4 PTFE eluent filter with 1 m PTFE tubing (ID 1.5 mm; OD 1/8")
  - 4 PTFE tubing pieces (ID 1.5 mm; OD 1/8") of 270, 280, 290, and 300 mm length
  - Connecting piece with 4 blind fittings
  - 30 Bushings
  - 30 Tefzel sealings
  - 1 Syringe with Luer-lock needle 1.5 x 50
- PEEK capillary AD1/16"x ID 0.25 mm
- 3 m Tefzel tubing AD1/16"x ID 0.7mm
- 15 Dynaseal bushings 1/16"
- 15 PETP double cone sealings
- Tube cutter
- Set of tools including
  - 1 wrench 1/4" x 5/16"
  - 1 hexagonal spanner 3 mm

## Technical Data

### Smartline Pump 1000

Delivery system	Double-piston pump with working and auxiliary piston
Flow rate range	0.001 – 9.999 ml/min
Flow rate accuracy	< 0.5 %, at 1 ml/min 12 MPa
Flow rate precision	< 0.1 %, at 1 ml/min 12 MPa
Residual pulsation	< 0.5 % at 1 ml/min methanol:water 8:2, 12 MPa
Max.operating pressure	40 MPa
Power supply	90 – 260 V, 47 – 63 Hz, 100 W
Dimensions	226 x 135 x 390 mm (W x H x D)
Weight	5.3 kg

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## Warranty statement

The warranty period of the Smartline Sample Preparation Unit 6300 is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Defective Smartline Sample Preparation Unit 6300 should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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## Declaration of conformity

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**Manufacturer's name and address:  
Herstellername und -adresse**

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**Smartline Sample Preparation Unit 6300,**

Order number, Bestellnummer **A5003**

Including / einschließlich

**Smartline Assistant 6010, V10**

complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)  
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)  
EN 6100-3-2 (10/2006)  
EN 61326-1 (10/2006)

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Dr. Ing. Herbert KNAUER GmbH